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CONTENTS

I. STATE SCIENTIFIC INSTITUTION “INSTITUTE OF PHYSIOLOGY OF THE NATIONAL ACADEMY OF SCIENCES OF BELARUS”	56
1. MULTILAYER VISUAL PLATFORM FOR ASSESSING CHANGES IN WOUND EXUDATE	56
II. STATE SCIENTIFIC INSTITUTION “INSTITUTE OF APPLIED PHYSICS OF THE NATIONAL ACADEMY OF SCIENCES OF BELARUS”, STATE SCIENTIFIC INSTITUTION “INSTITUTE OF PHYSIOLOGY OF THE NATIONAL ACADEMY OF SCIENCES OF BELARUS”	58
2. “BELCARDIOAI” — AN INTELLIGENT SYSTEM FOR ANALYZING DATA FROM AN EXPERIMENTAL CARDIAC RECORDER	58
III. STATE SCIENTIFIC INSTITUTION “B. I. STEPANOV INSTITUTE OF PHYSICS OF THE NATIONAL ACADEMY OF SCIENCES OF BELARUS”	60
3. AIR PLASMA JET GENERATOR “ALOE”	60
IV. STATE SCIENTIFIC INSTITUTION “INSTITUTE OF BIOORGANIC CHEMISTRY OF THE NATIONAL ACADEMY OF SCIENCES OF BELARUS”	61
4. BIOSPECIFIC HEMOSORBENT “ANTICYTOKINE-HEMO”	61
V. STATE SCIENTIFIC INSTITUTION “INSTITUTE OF MICROBIOLOGY OF THE NATIONAL ACADEMY OF SCIENCES OF BELARUS”	63
5. ENZYME PREPARATION LACTATE OXIDASE	63
VI. REPUBLICAN RESEARCH UNITARY ENTERPRISE “INSTITUTE OF BIOCHEMISTRY OF BIOLOGICALLY ACTIVE COMPOUNDS OF THE NATIONAL ACADEMY OF SCIENCES OF BELARUS”	64
6. BIOLOGICALLY ACTIVE FOOD SUPPLEMENT “IMMUNON-BIO”	64
7. TWO-COMPONENT POLYURETHANE COMPOUND FOR 3D PRINTING	64
8. MEDICAL DEVICE “RANLEK”	65
9. SOLUTION FOR ELIMINATION AND DECONTAMINATION OF DNA/RNA “BIOS”	66
10. HAIR GROWTH STIMULATION PRODUCT “CRINIS”	67
VII. JOINT-STOCK COMPANY “NPO CENTER” OF THE NATIONAL ACADEMY OF SCIENCES OF BELARUS, EDUCATIONAL INSTITUTION “THE BELARUSIAN STATE MEDICAL UNIVERSITY”	68
11. THE DEVICE FOR HYPERGRAVITATIONAL THERAPY IN THE CRANIO-CAUDAL DIRECTION GU-20.001.00.00.000	68
VIII. RESEARCH INSTITUTE OF HYGIENE, TOXICOLOGY, EPIDEMIOLOGY, VIROLOGY AND MICROBIOLOGY OF THE STATE INSTITUTION “REPUBLICAN CENTER OF HYGIENE, EPIDEMIOLOGY AND PUBLIC HEALTH”	69
12. METHOD FOR ASSESSING THE INTEGRAL TOXICITY OF AQUEOUS HEAVY METAL SOLUTIONS IN AN AGGRAVATED TEST USING THE TEST ORGANISM <i>CYPRIDOPSIS VIDUA</i>	69
13. METHOD FOR MEASURING THE MASS CONCENTRATION OF COBALT, NICKEL, ZINC, CADMIUM, IRON, MANGANESE, MOLYBDENUM, TIN, TUNGSTEN, VANADIUM, LEAD, COPPER, AND CHROMIUM IN WORKPLACE AIR USING INDUCTIVELY COUPLED PLASMA ATOMIC EMISSION SPECTROMETRY	70
14. METHOD OF MEASURING MASS CONCENTRATIONS OF ACRYLIC ACID NITRILE, METHYL ACRYLATE, AND METHYL METHACRYLATE IN THE AIR OF THE WORKING AREA	71
15. HEALTH RISK ASSESSMENT METHOD ASSOCIATED WITH MIGRATION OF CHEMICALS FROM MATERIALS INTENDED FOR CONTACT WITH FOOD PRODUCTS	72
16. METHOD OF HYGIENIC STANDARDIZATION OF THE CONTENT OF ANTIBACTERIAL DRUGS IN THE AIR OF THE WORKING AREA	73
17. METHODS OF HYGIENIC ASSESSMENT OF THE CHEMICAL COMPOSITION OF REPRESENTATIVE DIETS	74
IX. THE STATE INSTITUTION “REPUBLICAN SCIENTIFIC AND PRACTICAL CENTER FOR CHILDREN’S ONCOLOGY, HEMATOLOGY AND IMMUNOLOGY”	75
18. DEVELOPMENT OF A CAR RECEPTOR STRUCTURE FOR NK CELL MODIFICATION	75
19. BIOINFORMATICS ALGORITHM FOR AUTOMATED DETECTION AND INTERPRETATION OF GENETIC VARIANTS ASSOCIATED WITH PRIMARY IMMUNODEFICIENCIES	76

X. STATE INSTITUTION “REPUBLICAN SCIENTIFIC AND PRACTICAL CENTER FOR PULMONOLOGY AND TUBERCULOSIS”	77
20. METHOD FOR DETERMINING THE DRUG SENSITIVITY OF <i>MYCOBACTERIUM ABSCESSUS</i> TO MACROLIDES AND AMINOGLYCOSIDES	77
21. ALGORITHM FOR SELECTING A CHEMOTHERAPY REGIMEN FOR PATIENTS WITH MULTIDRUG-RESISTANT TUBERCULOSIS, INCLUDING FLUOROQUINOLONE RESISTANCE.....	78
22. INSTRUCTION DESCRIBING THE “METHODS FOR IDENTIFICATION OF THE <i>MYCOBACTERIUM TUBERCULOSIS</i> GENOTYPE”	79
XI. REPUBLICAN RESEARCH AND PRACTICAL CENTER OF NEUROLOGY AND NEUROSURGERY, REPUBLICAN INNOVATIVE UNITARY ENTERPRISE “SCIENCE AND TECHNOLOGY PARK OF BNTU “POLYTECHNIC”	81
23. MEDICAL DEVICE “SYNTHETIC PROTECTOR FOR MICROVASCULAR DECOMPRESSION”	81
XII. STATE INSTITUTION “REPUBLICAN SCIENTIFIC AND PRACTICAL CENTER FOR TRANSFUSIOLOGY AND MEDICAL BIOTECHNOLOGY”	82
24. THE DEVELOPMENT OF THE DOMESTIC MEDICINAL PRODUCT “PLAZMABEL” BASED ON VIRUS INACTIVATED CRYOPLASMA.....	82
XIII. THE EDUCATIONAL, SCIENTIFIC, AND PRODUCTION REPUBLICAN UNITARY ENTERPRISE “UNITECHPROM BSU”, THE BSU INSTITUTION “RESEARCH INSTITUTE OF PHYSICO-CHEMICAL PROBLEMS”	83
25. ORIGINAL ANTITUMOR DRUG “PROSPIDELONG”	83
XIV. EDUCATIONAL INSTITUTION “BELARUSIAN STATE MEDICAL UNIVERSITY”	84
26. SYNERGISTIC COMBINATION OF BEMITYL AND THYMOQUINONE: NEW HORIZONS IN THE FIGHT AGAINST HYPOXIA AND EXTREME CONDITIONS	84
XV. EDUCATIONAL INSTITUTION “GOMEL STATE MEDICAL UNIVERSITY”	86
27. ANTIBACTERIAL TUBULAR VASCULAR IMPLANT.....	86
XVI. VITEBSK STATE ORDER OF PEOPLES’ FRIENDSHIP MEDICAL UNIVERSITY	87
28. COMPREHENSIVE DIAGNOSIS OF METAL HYPERSENSITIVITY	87
29. DUAL-LUMEN NEPHROSTOMY CATHETER	88
30. REVERSE RETINACULOTOMY	89
31. ANDROID MOBILE APP “LIVERSCORE”	89
XVII. EDUCATIONAL INSTITUTION “BELARUSIAN STATE UNIVERSITY OF INFORMATICS AND RADIOELECTRONICS”, STATE INSTITUTION “REPUBLICAN SCIENTIFIC AND PRACTICAL CENTER FOR PULMONOLOGY AND PHTHISIOLOGY”	92
32. SYSTEM FOR THE AUTOMATED DETERMINATION AND PERSONALIZED ADJUSTMENT OF OXYGEN THERAPY PARAMETERS IN REAL TIME BASED ON THE PATIENT’S CURRENT CONDITION	92
XVIII. APEX-REAL LIMITED LIABILITY COMPANY (RESIDENT OF BREST SCIENCE AND TECHNOLOGY PARK)	95
33. DEVELOPMENT AND IMPLEMENTATION OF A NEW ELECTRONIC SYSTEM FOR MANAGING AND DISTRIBUTING PATIENT FLOWS	95
XIX. FEDERAL STATE BUDGETARY EDUCATIONAL INSTITUTION OF HIGHER EDUCATION “KEMEROVO STATE MEDICAL UNIVERSITY” OF THE MINISTRY OF HEALTH OF THE RUSSIAN FEDERATION	96
34. A DEVICE FOR SELF-COLLECTION OF VAGINAL DISCHARGE — THE “JUNO” PROBE	96

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E-mail: artymuk@gmail.com

MEDICINE AND PHARMACY

I. STATE SCIENTIFIC INSTITUTION “INSTITUTE OF PHYSIOLOGY OF THE NATIONAL ACADEMY OF SCIENCES OF BELARUS”

1. MULTILAYER VISUAL PLATFORM FOR ASSESSING CHANGES IN WOUND EXUDATE

BRIEF DESCRIPTION OF INNOVATIVE DEVELOPMENT

The innovative development is a multilayer chitosan-containing visual test system for semiquantitative assessment of local changes in wound exudate. It is intended for supplementary monitoring of the wound environment during dressing changes without laboratory equipment, electronics, a power source or a dedicated reader. The system is a disposable multilayer patch comprising an atraumatic contact layer, protective separation membranes, a chitosan-containing absorption and distribution layer, sensor zones and a transparent outer covering. Exudate reaches the sensor elements and produces visually detectable colour changes. The basic prototype includes a pH sensor, a protease sensor and a filling control zone for assessing changes in acid-base conditions, proteolytic activity and exudate inflow. Its autonomous operation, simple visual readout and limited infrastructure requirements make the development promising for field and mobile medical settings, first aid and premedical care, remote medical posts, outpatient practice and home wound care. A laboratory demonstration prototype has been produced and tested using model solutions.

A clean prototype before the reaction/contact

Hemostasis/remodeling

Inflammation/exudation

Proliferation/granulation

TECHNICAL ADVANTAGES

The development belongs to the current field of smart wound dressings, which convert wound environment parameters into visual, optical or electrical signals. International solutions in this class are predominantly focused on electronic and wearable sensor systems requiring readers, power sources, protection of electronic components and stable signals. The proposed development combines assessment of multiple wound exudate parameters with an entirely passive visual operating principle. The basic prototype assesses pH, protease activity and exudate filling without electronics or specialised equipment. Results can be interpreted directly from changes in the colour of the sensor zones. An additional design advantage is the multilayer architecture, incorporating a chitosan-containing absorption and distribution element and barriers separating sensor reagents from the wound surface, with directed exudate transport towards the diagnostic zones. This approach could reduce the cost and technical complexity of monitoring and is particularly promising for field, mobile and resource-limited settings.

EXPECTED RESULT OF APPLICATION

The development is expected to enable rapid visual semiquantitative assessment of changes in the wound environment during dressing changes without laboratory equipment, electronics or power sources. The system is intended to track changes in pH, proteolytic activity and exudation intensity, potentially supporting earlier detection of adverse changes and informed adjustment of wound management. Its greatest practical value is anticipated in field, mobile and resource-limited settings where access to laboratory testing is restricted.

CURRENT STAGE OF DEVELOPMENT

A laboratory demonstration prototype of the multilayer visual test system has been produced. The design principle, model exudate transport and visual responses of the sensor zones have been tested using model solutions. Clinical evaluation and further validation have not been performed.

INFORMATION ON THE LEGAL PROTECTION

Intellectual property protection has not yet been formalised. Materials are being prepared to determine the appropriate form of protection and subsequently seek patent protection for the technical solution.

POTENTIAL CONSUMERS AND INTERESTED PARTIES

Medical organisations in the Republic of Belarus and abroad: outreach and mobile medical services, field medical units, remote medical posts, surgical and outpatient facilities, and community nursing and wound care services.

CONTACT INFORMATION

Palina Litvinava, Researcher and Postgraduate Researcher.

Phone number: (+375 29) 185 28 88

E-mail: palinalitvinava@gmail.com

Maksim Simirski, Acting Head of General Practice Department No. 2 at the 19th Central District Polyclinic, Postgraduate Researcher.


Phone number: (+375 29) 158 71 77

E-mail: simirski98@gmail.com

II. STATE SCIENTIFIC INSTITUTION “INSTITUTE OF APPLIED PHYSICS OF THE NATIONAL ACADEMY OF SCIENCES OF BELARUS”, STATE SCIENTIFIC INSTITUTION “INSTITUTE OF PHYSIOLOGY OF THE NATIONAL ACADEMY OF SCIENCES OF BELARUS”

2. “BELCARDIOAI” — AN INTELLIGENT SYSTEM FOR ANALYZING DATA FROM AN EXPERIMENTAL CARDIAC RECORDER

BRIEF DESCRIPTION OF INNOVATIVE DEVELOPMENT

Field of application: scientific and clinical-physiological research on the cardiovascular system.

Purpose: “BelCardioAI” intelligent system is designed to perform intelligent analysis of recordings comprising three electrocardiogram (ECG) channels and one photoplethysmogram (PPG) channel, and to generate a clinical conclusion and a PDF report.

Key design, technological, and operational characteristics: the experimental cardiac recorder for which this project was developed is an original domestic product lacking certified commercial software for ECG and PPG analysis; key device parameters (such as sampling rate and lead orientation) were undocumented by the manufacturer and had to be determined and validated during the project. “BelCardioAI” addresses this limitation by providing a complete software pipeline ranging from reading raw analog-to-digital converter samples to generating a formalized clinical conclusion in Russian. The device used in the project is an early hardware version; a portable variant of the system already records PPG, blood oxygen saturation, and body temperature, while later hardware versions are slated to include bioimpedance and body position monitoring. Extending the AI pipeline to incorporate these additional parameters of the full system is an anticipated outcome of the project’s next stage, rather than a currently implemented feature.

TECHNICAL ADVANTAGES

2 trained AI components: a classifier based on a 1D convolutional neural network for categorizing cardiac complexes, and an unsupervised autoencoder trained directly on the recording being analyzed to detect atypical segments without prior labeling; both components operate entirely locally, without network access.

An AI integrity principle: any model failing to outperform a trivial baseline prediction is automatically disabled and its output suppressed — for instance, blood pressure estimation is automatically disabled for the current dataset if it does not meet this standard.

Classical digital signal processing algorithms serving as the foundation for the AI layer: channel filtering and quality assessment, an adaptive Pan-Tompkins algorithm for R-wave detection (with multi-channel gap recovery), automatic annotation of P-waves, QRS complexes, and T-waves, and lead polarity normalization.

Calculation of over 40 heart rate variability (HRV) metrics — including time-domain, frequency-domain, and non-linear parameters, as well as R. M. Baevsky's variational pulsometry indices.

Automatic classification of rhythm and conduction disorders (extrasystoles, pauses, signs of atrial fibrillation, AV-block, QTc prolongation) with cross-verification of findings across multiple leads simultaneously.

Automatic generation of structured clinical findings and PDF reports in Russian, along with cohort analysis across the entire study population.

A fully autonomous (offline) desktop version — a Windows executable requiring no internet connection, Python installation, or reliance on cloud-based AI services or API keys, a critical requirement for handling medical data.

EXPECTED RESULT OF APPLICATION

Reduction in the time required to analyze a single cardiac recording (taking just a few seconds per session).

Standardization and reproducibility of analysis results for an experimental device that previously lacked dedicated, expert-grade software.

Expanded diagnostic and research capabilities through automated screening for rhythm and conduction disorders, heart rate variability, signs of atrial fibrillation, and apnea episodes (detected via ECG-derived respiration).

A technical foundation for the future certification of the domestic cardiac recorder and its software as a medical device.

Technological independence from foreign cloud-based AI services — all AI components operate locally without transmitting data outside the organization, a crucial factor for import substitution and medical data security requirements.

CURRENT STAGE OF DEVELOPMENT

A research project currently in further development, specifically regarding the creation of software for the real-time automated analysis of incoming data.

INFORMATION ON THE LEGAL PROTECTION

No.

POTENTIAL CONSUMERS AND INTERESTED PARTIES

Healthcare and sports medicine institutions — state (including departmental) and private.

CONTACT INFORMATION

Petr Khmarsky, Leading Researcher at the Laboratory of Microwave Electronics and Radiotomography, Candidate of Technical Sciences, Associate Professor.

Phone number: (+375 29) 219 56 49

E-mail: pierre2009@mail.ru

Yegor Lemesenko, Head of the Center of Aerospace Biomedicine and Rehabilitation, Candidate of Medical Sciences, Associate Professor.

Phone number: (+375 29) 621 12 71

E-mail: iflemeshko@gmail.com

III. STATE SCIENTIFIC INSTITUTION “B. I. STEPANOV INSTITUTE OF PHYSICS OF THE NATIONAL ACADEMY OF SCIENCES OF BELARUS”

3. AIR PLASMA JET GENERATOR “ALOE”

BRIEF DESCRIPTION OF INNOVATIVE DEVELOPMENT

“ALOE” is a device that generates an air plasma jet from atmospheric air, containing biologically active gases (NO, NO₂ and HNO₂) with high antimicrobial activity. Fields of application — medicine and veterinary medicine.

TECHNICAL ADVANTAGES

The effectiveness of “ALOE” device is ensured by the generation of biologically active gases (NO, NO₂ and HNO₂) from atmospheric air, which possess high antimicrobial properties. The device operates from a 220 V / 50 Hz mains supply, uses atmospheric air, features electronic automated control, is mobile, and weighs 8 kg. In terms of resistance to mechanical impacts and electrical safety, the device meets the requirements for electronic devices for medical and laboratory use. In terms of resistance to climatic impacts, it complies with conditions with temperature limits from –10 to +35 °C and humidity up to 80 %. The degree of protection corresponds to IP54 according to IEC 60529-2013. It has no domestic analogues.

EXPECTED RESULT OF APPLICATION

The device is recommended for use in the complex therapy of dental diseases to provide antimicrobial, anti-inflammatory and regenerative effects, as well as for use in medical practice in healthcare institutions of the Republic of Belarus. In terms of its technical parameters, functional purpose and design, it meets the requirements imposed by medical science and clinical practice on physiotherapeutic equipment.

CURRENT STAGE OF DEVELOPMENT

10 prototype samples have been produced.

INFORMATION ON THE LEGAL PROTECTION

Method for producing bioactive nitrogen components and device for its implementation: Pat. 24764 Rep. of Belarus: IPC H 05H 1/00, A 61M 11/00 / L. V. Simonchik, A. V. Kazak; applicant and patent holder: State Scientific Institution “B. I. Stepanov Institute of Physics of the National Academy of Sciences of Belarus”. — No. a20240222; filed 22.10.2024; published 20.12.2025.

POTENTIAL CONSUMERS AND INTERESTED PARTIES

Healthcare institutions of the Ministry of Health of the Republic of Belarus.

CONTACT INFORMATION

Leonid Simonchik, Chief Researcher, Doctor of Physical and Mathematical Sciences, Professor.

Phone number: (+375 29) 150 62 06

E-mail: l.simonchik@ifanbel.bas-net.by

Alexandra Kazak, Leading Researcher, Candidate of Physical and Mathematical Sciences, Associate Professor.

Phone number: (+375 29) 257 00 71

E-mail: a.pavlova@ifanbel.bas-net.by

IV. STATE SCIENTIFIC INSTITUTION “INSTITUTE OF BIOORGANIC CHEMISTRY OF THE NATIONAL ACADEMY OF SCIENCES OF BELARUS”

4. BIOSPECIFIC HEMOSORBENT “ANTICYTOKINE-HEMO”

BRIEF DESCRIPTION OF INNOVATIVE DEVELOPMENT

Excessive concentrations of pro-inflammatory cytokines in the blood can sustain and exacerbate systemic inflammation, contributing to the development of severe pathological conditions. Current therapies primarily focus on suppressing cytokine production or activity through pharmacological agents.

Our proposed approach involves directly removing excess target cytokines from the bloodstream rather than blocking the inflammatory signal.

The anti-cytokine hemosorbent under development consists of a biocompatible sorption material featuring functional ligands capable of selectively binding pro-inflammatory cytokines directly from the blood.

The technology relies on a combination of polymer materials engineering, peptide ligand design, and molecular modeling, utilizing state-of-the-art approaches. The ligands are engineered to ensure high affinity for the target cytokines while maintaining the essential properties of the hemocompatible sorption material.

TECHNICAL ADVANTAGES

Selectivity. Instead of the non-specific removal of a large number of blood components, the technology focuses on interacting with specific target molecules.

Direct action on the source of the pathological inflammatory signal. The hemosorbent acts directly on circulating cytokines, potentially enabling a rapid reduction in their concentration during the extracorporeal blood purification procedure.

Modular platform. Modifying the ligand structure potentially allows for the creation of sorbents targeting either individual cytokines or various combinations of targets.

Compatibility with hemosorption technology. The development aims to produce a sorbent material suitable for use in specialized mass-transfer devices for extracorporeal therapy.

Potential for a personalized approach. In the future, the technology could serve as a basis for selecting sorbents tailored to the specific profile of inflammatory mediators circulating in a patient.

Commercial potential. The project aims to create a new generation of selective hemosorbents for use in extracorporeal therapy for conditions characterized by pathologically elevated levels of pro-inflammatory cytokines.

The technology's platform-based design enables the creation of a product line based on a common polymer matrix and peptide ligands capable of binding different classes of cytokines.

This creates the potential to develop not just a single medical device, but a technological foundation for selective blood purification.

Potential design of the hemoadsorbent's appearance

EXPECTED RESULT OF APPLICATION

The successful application of this development is expected to result, at a minimum, in an expanded arsenal of treatments for conditions associated with elevated levels of pro-inflammatory cytokines in the blood and the active integration of sorbents into therapeutic protocols in the Republic of Belarus; at a maximum, it is expected to lead to market entry and a strong foothold in the EAEU, the Middle East, Africa, and other partner nations of the Republic of Belarus.

CURRENT STAGE OF DEVELOPMENT

An experimental sample has been obtained.

INFORMATION ON THE LEGAL PROTECTION

Trade secrets.

POTENTIAL CONSUMERS AND INTERESTED PARTIES

Organizations under the Ministry of Health of the Republic of Belarus; organizations under the ministries of health of purchasing countries.

CONTACT INFORMATION

Ivan Semkov, Junior Researcher at the Laboratory of Applied Biochemistry.

Phone number: (+375 25) 910 35 23

E-mail: semkov_i@mail.ru

V. STATE SCIENTIFIC INSTITUTION “INSTITUTE OF MICROBIOLOGY OF THE NATIONAL ACADEMY OF SCIENCES OF BELARUS”

5. ENZYME PREPARATION LACTATE OXIDASE

BRIEF DESCRIPTION OF INNOVATIVE DEVELOPMENT

The lactate oxidase enzyme catalyzes the oxidation of L-lactate to pyruvate, producing hydrogen peroxide; it can be used as a bioreagent in enzymatic biosensors and analytical kits for the quantitative determination of lactate in biological fluids. Its field of application is clinical diagnostics: assessment of tissue hypoxia, monitoring of critical conditions (sepsis, shock, respiratory failure), neonatology, and sports medicine. A key characteristic of the enzyme is its high substrate specificity for L-lactate; the measurement range covers clinically significant concentrations. Technical and operational indicators include high selectivity, result reproducibility.

TECHNICAL ADVANTAGES

The development meets global standards regarding analytical characteristics and surpasses foreign counterparts in terms of operational stability and non-invasive monitoring capabilities. The project addresses the need for import substitution of a critical bioreagent and establishes a foundation for producing competitive domestic diagnostic solutions.

EXPECTED RESULT OF APPLICATION

Adopting the domestic bioreagent will reduce reliance on imported enzyme preparations and consumables, while the scalable production technology — covering everything from producer cultivation to lyophilization and biochip fabrication — will enable the serial production of sensors and laboratory kits.

CURRENT STAGE OF DEVELOPMENT

Research or development (technological) work has been completed.

INFORMATION ON THE LEGAL PROTECTION

No.

POTENTIAL CONSUMERS AND INTERESTED PARTIES

Manufacturers of biosensor systems.

CONTACT INFORMATION

Semashko Tatiana, Head of the Micromycete Biotechnology Laboratory, Candidate of Biological Sciences, Associate Professor.

Phone number: (+375 17) 236 50 63

E-mail: tsemashko@mbio.bas-net.by

VI. REPUBLICAN RESEARCH UNITARY ENTERPRISE “INSTITUTE OF BIOCHEMISTRY OF BIOLOGICALLY ACTIVE COMPOUNDS OF THE NATIONAL ACADEMY OF SCIENCES OF BELARUS”

6. BIOLOGICALLY ACTIVE FOOD SUPPLEMENT “IMMUNON-BIO”

BRIEF DESCRIPTION OF INNOVATIVE DEVELOPMENT

“IMMUNON-BIO” is a biologically active food supplement with a complex positive effect on the body. The product helps strengthen the immune system, supports stable liver function and reduces the risk of side effects from the gastrointestinal tract.

TECHNICAL ADVANTAGES

The unique combination of components in “IMMUNON-BIO” ensures the intake of an essential trace element required for the formation of proteins and amino acids. This trace element plays a key role in the synthesis of endogenous enzymes and hormones that regulate metabolic processes, serves as one of the main nutrients of the immune system, and is a component of insulin, helping to prolong its hypoglycemic effect. In addition, it potentiates the action of testosterone and folliculin, stimulates the formation of red blood cells and hemoglobin, and also has lipotropic properties.

EXPECTED RESULT OF APPLICATION

With regular intake of “IMMUNON-BIO”, the following is expected: strengthening of immunity and reduction of the risk of colds and viral diseases; tissue regeneration; destruction of proteins produced by bacteria; maintenance of male reproductive system health; protection of pancreatic cells that produce insulin; stable liver function; improvement of the condition of skin, nails and hair (prevents acne and hair loss); synthesis of the body's own collagen; activation of the functions of the organs of vision, taste and smell.

CURRENT STAGE OF DEVELOPMENT

The development has been introduced into production.

INFORMATION ON THE LEGAL PROTECTION

Technical specifications, trademark, EAC certificate, registration certificate of the Ministry of Health, sanitary-hygienic registration.

POTENTIAL CONSUMERS AND INTERESTED PARTIES

Ministry of Health of the Republic of Belarus, the general public, sanatorium-resort institutions.

CONTACT INFORMATION

Oleg Kuznetsov, Director, Candidate of Biological Sciences, Associate Professor.
Phone number: (+375 15) 255 90 75
E-mail: office@ibiochemistry.by

7. TWO-COMPONENT POLYURETHANE COMPOUND FOR 3D PRINTING

BRIEF DESCRIPTION OF INNOVATIVE DEVELOPMENT

A two-component polyurethane (polymer) cold-curing compound intended for 3D printing. The material makes it possible to create three-dimensional models of objects in three-dimensional space (prototyping)

and to manufacture structural parts for various industries. A separate area of application is the production of medical products, including means of individual immobilization and fixation of body parts for traumatology, correction and rehabilitation.

TECHNICAL ADVANTAGES

The key advantage of the compound is cold curing, which does not require heating or special equipment, simplifying and reducing the cost of the production process. The two-component composition ensures precise control of the properties of the finished product and stability of characteristics, while the polyurethane base combines strength and elasticity, which is especially important for products in contact with the human body. The material is versatile: it is suitable both for rapid prototyping and for the manufacture of functional structural parts, including individual medical products tailored to the anatomical features of a specific patient.

EXPECTED RESULT OF APPLICATION

Creation of individual immobilizing and fixing products for various body parts, used in traumatology, correction and rehabilitation, with a high level of comfort and precision of fitting to the patient.

CURRENT STAGE OF DEVELOPMENT

The development has been introduced into production.

INFORMATION ON THE LEGAL PROTECTION

Technical specifications, trademark, EAC certificate, registration certificate of the Ministry of Health, sanitary-hygienic registration.

POTENTIAL CONSUMERS AND INTERESTED PARTIES

Ministry of Health of the Republic of Belarus, the general public, sanatorium-resort institutions.

CONTACT INFORMATION

Oleg Kuznetsov, Director, Candidate of Biological Sciences, Associate Professor.
Phone number: (+375 15) 255 90 75
E-mail: office@ibiochemistry.by

8. MEDICAL DEVICE "RANLEK"

BRIEF DESCRIPTION OF INNOVATIVE DEVELOPMENT

A medical device used as an auxiliary remedy for sunburn and thermal burns, as well as for skin recovery in various injuries: wounds, burns, trophic and radiation ulcers. The product is used for the prevention of bacterial complications and local anesthesia, and is also suitable for daily skin care, including very dry and flaking-prone skin.

TECHNICAL ADVANTAGES

"RANLEK" is the only composition on the EAEU market that simultaneously combines wound-healing, analgesic and antibacterial effects.

EXPECTED RESULT OF APPLICATION

When used, "RANLEK" improves skin condition, provides protection against the negative effects of the environment, has a soothing and softening effect, and also helps reduce redness. Due to its properties, the product is suitable for first-aid kits.

CURRENT STAGE OF DEVELOPMENT

The development has been introduced into production.

Medical device "RANLEK"

INFORMATION ON THE LEGAL PROTECTION

Technical specifications, trademark, EAC certificate, registration certificate of the Ministry of Health, sanitary-hygienic registration.

POTENTIAL CONSUMERS AND INTERESTED PARTIES

Ministry of Health, healthcare institutions (hospitals, polyclinics, trauma and burn centers), pharmacy chains, sanatorium-resort organizations, emergency medical services, as well as the general public for home use and daily skin care.

CONTACT INFORMATION

Oleg Kuznetsov, Director, Candidate of Biological Sciences, Associate Professor.

Phone number: (+375 15) 255 90 75

E-mail: office@ibiochemistry.by

9. SOLUTION FOR ELIMINATION AND DECONTAMINATION OF DNA/RNA "BIOS"

BRIEF DESCRIPTION OF INNOVATIVE DEVELOPMENT

"BIOS" is a solution designed for the elimination of RNases and DNases, as well as the decontamination of DNA from various surfaces. The product ensures effective cleaning of work areas and equipment from traces of nucleic acids and enzymes, which is critically important for obtaining reliable results in molecular biological research.

TECHNICAL ADVANTAGES

"BIOS" is the first domestic product of this purpose. Its unique composition allows for the treatment of PCR equipment, electrophoresis equipment, as well as consumables. The product is completely ready for use and requires no additional preparation.

EXPECTED RESULT OF APPLICATION

The use of "BIOS" provides an economic effect 2.8–5.2 times greater compared to any imported analogue, making it a cost-effective solution for laboratories. The product is ready for use, which reduces labor costs and time for preparing working solutions.

CURRENT STAGE OF DEVELOPMENT

The development has been introduced into production.

INFORMATION ON THE LEGAL PROTECTION

Technical specifications, trademark, EAC certificate, registration certificate of the Ministry of Health, sanitary-hygienic registration.

POTENTIAL CONSUMERS AND INTERESTED PARTIES

Ministry of Health, research institutes and laboratories of molecular biological profile, PCR laboratories, centers of hygiene and epidemiology, veterinary and food laboratories, educational institutions (departments of biology, medicine, genetics), as well as manufacturers of diagnostic kits and consumables for molecular biology.

CONTACT INFORMATION

Oleg Kuznetsov, Director, Candidate of Biological Sciences, Associate Professor.

Phone number: (+375 15) 255 90 75

E-mail: office@ibiochemistry.by

10. HAIR GROWTH STIMULATION PRODUCT “CRINIS”

BRIEF DESCRIPTION OF INNOVATIVE DEVELOPMENT

“Crisis” is a 100 % natural product in the form of a spray for external use. Its active base is a water-alcohol extract obtained from plant raw materials. The product is intended for both men and women suffering from increased hair loss. Its mechanism of action is aimed at strengthening hair follicles, preventing further hair loss and actively stimulating the growth of new hair.

TECHNICAL ADVANTAGES

The key advantage of the development is its natural composition, which eliminates the risk of systemic side effects characteristic of many synthetic analogues. The product effectively addresses two main objectives: it stops the process of hair loss and triggers the growth of new hair, acting directly on the root system.

EXPECTED RESULT OF APPLICATION

With regular use over the course of a month, the product reduces the intensity of hair loss by up to 85 %. There is an improvement in the health of the hair bulb (it becomes shiny and viable), as well as restoration of the shine and structure of hair along its entire length — from roots to ends.

CURRENT STAGE OF DEVELOPMENT

The development has been introduced into production.

INFORMATION ON THE LEGAL PROTECTION

Technical specifications, trademark, EAC certificate, registration certificate of the Ministry of Health, sanitary-hygienic registration.

POTENTIAL CONSUMERS AND INTERESTED PARTIES

Ministry of Health of the Republic of Belarus, the general public, sanatorium-resort institutions.

CONTACT INFORMATION

Oleg Kuznetsov, Director, Candidate of Biological Sciences, Associate Professor.
Phone number: (+375 15) 255 90 75
E-mail: office@ibiochemistry.by

Product “Crisis”

VII. JOINT-STOCK COMPANY “NPO CENTER” OF THE NATIONAL ACADEMY OF SCIENCES OF BELARUS, EDUCATIONAL INSTITUTION “THE BELARUSIAN STATE MEDICAL UNIVERSITY”

11. THE DEVICE FOR HYPERGRAVITATIONAL THERAPY IN THE CRANIO-CAUDAL DIRECTION GU-20.001.00.00.000

BRIEF DESCRIPTION OF INNOVATIVE DEVELOPMENT

The unit is designed for hypergravitational physiotherapy in the cranio-caudal direction. The method is based on the beneficial effect of gravity (G) on the human body. The principle of operation consists in applying an overload (g) to the patient's body, positioned horizontally, through the rotation of a platform. Centrifugal force creates a vector of longitudinal overload (gy) directed from head to feet, increasing the force of gravity (Gx) in the prone direction (from back to chest). This causes a complex of positive physiological reactions: improvement of regional blood flow and microcirculation, restoration of transcapillary exchange, acceleration of tissue regeneration, and elimination of structural deformations in joints. The unit has undergone a full cycle of testing: laboratory and technical examinations, acceptance and clinical trials at the Healthcare Institution “4th City Clinical Hospital named after N. E. Savchenko” and the Educational Institution “BSMU” (March 2025). A conclusion was obtained from the Republican Unitary Enterprise “Center for Expertise and Testing in Healthcare” (April 2025) and an indefinite registration certificate No. IM-7.117169 (May 2025).

Main technical characteristics: rated voltage $230\text{ V} \pm 23\text{ V}$, $50\text{ Hz} \pm 0.5\text{ Hz}$; maximum power consumption 1.8 kW; dimensions $2700 \times 1300 \times 1500\text{ mm}$; weight up to 300 kg; box load capacity up to 145 kg; patient height from 140 to 190 cm; drive performance up to 47 rpm; dosed exposure to simulated artificial gravity in the cranio-caudal direction up to 3.5 g.

TECHNICAL ADVANTAGES

The unique development is an adaptive technology for recovery and treatment using the method of hypergravitational physiotherapy. It is based on the organization of the beneficial effect of gravity in the cranio-caudal direction, which distinguishes this method from traditional physiotherapeutic approaches. The unit allows for a complex effect on the body, ensuring improved microcirculation, accelerated regeneration and restoration of transcapillary exchange, which has been confirmed by clinical trials. No direct analogues implementing this principle in such a design have been identified on the market.

EXPECTED RESULT OF APPLICATION

The introduction of the unit into medical practice will improve the quality of treatment and rehabilitation of patients with a wide range of diseases: in traumatology and orthopedics, surgery, cardiology, urology, gynecology and otorhinolaryngology. An increase in the effectiveness of therapy is expected due to non-invasive exposure, a reduction in recovery time and an improvement in the general physiological state of patients.

CURRENT STAGE OF DEVELOPMENT

A prototype sample has been produced.

INFORMATION ON THE LEGAL PROTECTION

No.

POTENTIAL CONSUMERS AND INTERESTED PARTIES

Healthcare institutions of the Republic of Belarus (hospitals, clinics, rehabilitation centers), as well as sanatorium-resort institutions.

CONTACT INFORMATION

Igor Kovalevsky, Director of the Branch “Center-5”.

Phone number: (+375 17) 259 03 15

E-mail: ikavaleuski@npo-center.by

VIII. RESEARCH INSTITUTE OF HYGIENE, TOXICOLOGY, EPIDEMIOLOGY, VIROLOGY AND MICROBIOLOGY OF THE STATE INSTITUTION “REPUBLICAN CENTER OF HYGIENE, EPIDEMIOLOGY AND PUBLIC HEALTH”

12. METHOD FOR ASSESSING THE INTEGRAL TOXICITY OF AQUEOUS HEAVY METAL SOLUTIONS IN AN AGGRAVATED TEST USING THE TEST ORGANISM *CYPRIDOPSIS VIDUA*

BRIEF DESCRIPTION OF INNOVATIVE DEVELOPMENT

Scope: preventive medicine, environmental hygiene, ecotoxicology.

Purpose: the method can be used within a suite of medical services aimed at the medical prevention of the adverse effects of toxicants, to establish a battery of test models based on the analysis of results from the quantitative assessment of the integral toxicity of aqueous media.

Main characteristics: express method for the quantitative assessment of the integral toxicity of solutions containing toxicants — including heavy metals — using an intensified test that involves exposing the crustacean *Cypridopsis vidua*, recording the effect in terms of immobilization, and subsequently calculating an integral toxicity index (IT), with the testing conducted within a temperature range of 26–28 °C over a period of 48 hours.

TECHNICAL ADVANTAGES

Technical advantages: Criteria-based quantitative parameters have been established for assessing the adverse effects of toxicant solutions based on the integral toxicity (IT) index: an IT value of less than 25 % is classified as low integral toxicity; 25 to 50 % as moderate; 50 to 75 % as pronounced; and over 75 % as severe.

Scientific and technical level: corresponds to the VI technological order.

EXPECTED RESULT OF APPLICATION

The method expands the range of sensitive techniques for assessing the adverse effects of environmental toxicants and pollutants and reduces the determination time from 96 to 48 hours, justifying its recommendation as a rapid method.

CURRENT STAGE OF DEVELOPMENT

At the implementation stage.

INFORMATION ON THE LEGAL PROTECTION

Patent No. BY 24650.

POTENTIAL CONSUMERS AND INTERESTED PARTIES

Institutions carrying out state sanitary supervision, health care organizations, state medical and scientific organizations, state educational institutions carrying out training, advanced training and (or) retraining of specialists with higher or specialized secondary medical education in the Republic of Belarus and abroad.

CONTACT INFORMATION

Natalya Dudchik, Chief Researcher at the Laboratory for the Study of Microbiota in Human Environmental Objects and Molecular Biological Research, Doctor of Biological Sciences, Professor.

Phone number: (+375 17) 379 13 85

E-mail: micro_sanitary@rspch.by

13. METHOD FOR MEASURING THE MASS CONCENTRATION OF COBALT, NICKEL, ZINC, CADMIUM, IRON, MANGANESE, MOLYBDENUM, TIN, TUNGSTEN, VANADIUM, LEAD, COPPER, AND CHROMIUM IN WORKPLACE AIR USING INDUCTIVELY COUPLED PLASMA ATOMIC EMISSION SPECTROMETRY

The high-tech equipment used

BRIEF DESCRIPTION OF INNOVATIVE DEVELOPMENT

Scope of application: for use by laboratories conducting sanitary monitoring of air safety in the working zone of industrial enterprises and other organization.

Use: quantitative determination of the mass concentration of 13 toxic metals in the workplace air.

Main characteristics: the measurement technique is based on collecting airborne particulate of workplace air on the AFA-ChA-20 filters (manufactured by Kimry Gorky Factory JSC, Russia), subjecting the filters to "wet" acid mineralization to bring the target metals into solution, and measuring the atomic emission generated when the sample is nebulized into an argon plasma inductively excited by a radio-frequency electromagnetic field. The emission of metal atoms is directly proportional to the mass concentration of these elements in the test solution. The mass concentration of the metals is calculated from the emission measurement results using a linear calibration curve constructed with five calibration solutions for each metal.

TECHNICAL ADVANTAGES

Technical advantages: rapid air sampling in the workplace (10 min); use of non-explosive gases for measurements (argon instead of acetylene/propane); favorable metrological characteristics (for all determined elements, the repeatability limit ranges from 10 to 26 %, the intermediate precision limit from 12 to 26 %, and the relative expanded uncertainty from 11 to 25 %); simultaneous determination of 13 hygienic significance metals in workplace air.

Scientific and technical level: exceeds that of the best domestic analogues.

EXPECTED RESULT OF APPLICATION

The application of this measurement technique in state sanitary supervision practice will contribute to enhancing occupational safety at industrial enterprises and organizations by yielding reliable results for the measurement of the mass concentration of 13 toxic metals in workplace air.

CURRENT STAGE OF DEVELOPMENT

At the implementation stage.

INFORMATION ON THE LEGAL PROTECTION

No.

POTENTIAL CONSUMERS AND INTERESTED PARTIES

Institutions carrying out state sanitary supervision, health care organizations, state medical and scientific organizations, state educational institutions carrying out training, advanced training and (or) retraining of specialists with higher or specialized secondary medical education in the Republic of Belarus and abroad.

CONTACT INFORMATION

Anna Kuzovkova, Head of the Spectrometric Research Laboratory, PhD.

Phone number: (+375 17) 373 95 00

E-mail: zav_lsi@rspch.by

14. METHOD OF MEASURING MASS CONCENTRATIONS OF ACRYLIC ACID NITRILE, METHYL ACRYLATE, AND METHYL METHACRYLATE IN THE AIR OF THE WORKING AREA

BRIEF DESCRIPTION OF INNOVATIVE DEVELOPMENT

Scope of application: air control of the working area.

Purpose: quantitative determination of acrylic acid nitrile, methyl acrylate and methyl methacrylate in the air of the work area.

Main characteristics: the principle of the method is based on the concentration of acrylic acid nitrile, methyl acrylate and methyl methacrylate from the air of the working area onto sorption tubes filled with a solid sorbent, their subsequent thermal desorption, then gas chromatographic separation of the substances to be determined on two parallel quartz capillary columns with a fixed liquid phase deposited on them, their identification by retention times on two channels of PID and quantitative determination by the method of absolute calibration.

A metrologically certified measurement technique has been developed, which is characterized by selectivity and high sensitivity (the lower limit of quantitative determination for each substance is 0.1 mg/m^3), the presence of acceptable metrological indicators (repeatability limit is $< 20 \%$, intermediate precision limit is $< 30 \%$ and relative extended uncertainty is $< 35 \%$) for all substances.

TECHNICAL ADVANTAGES

Technical advantages: the use of a metrologically certified methodology will improve the quality and efficiency of monitoring the content of acrylic acid nitrile, methyl acrylate and methyl methacrylate in the air of the work area, and use widespread equipment equipped with testing laboratories that monitor workplaces.

Scientific and technical level: exceeds the level of the best domestic analogues.

EXPECTED RESULT OF APPLICATION

The use of a metrologically certified methodology will prevent damage to human health and minimize the negative impact on human health by increasing the degree of sanitary reliability of the air in the working area at industrial enterprises, which will contribute to reducing the risks of diseases caused by air quality and lead to an increase in the standard of living of the population.

CURRENT STAGE OF DEVELOPMENT

Research or development (technological) work has been completed.

INFORMATION ON THE LEGAL PROTECTION

No.

POTENTIAL CONSUMERS AND INTERESTED PARTIES

The measurement methodology will be used by bodies and institutions that carry out state sanitary supervision and occupational safety control, state medical scientific health organizations and other organizations.

CONTACT INFORMATION

Tatiana Krymskaya, Head of the Chromatographic Research Laboratory.

Phone number: (+375 17) 379 08 57

E-mail: chromatographic@rspch.by

15. HEALTH RISK ASSESSMENT METHOD ASSOCIATED WITH MIGRATION OF CHEMICALS FROM MATERIALS INTENDED FOR CONTACT WITH FOOD PRODUCTS

BRIEF DESCRIPTION OF INNOVATIVE DEVELOPMENT

Scope: the method can be used in a complex of medical services aimed at medical prevention of negative impacts of substances migrating from materials in contact with food products on consumer's health.

Purpose: the method allows to ensure the application of risk-oriented approach in justification and updating of hygienic safety indicators of materials (including new types) and standards of chemical substances, justification of management decisions in the field of hygienic safety of materials (in terms of the need to establish or revise hygienic standards, laboratory monitoring of materials safety, determining the area of use of materials in certain areas of food production, etc.); introduction of measures to ensure hygienic safety of materials (including new types of materials) and chemical substances.

Main characteristics: the method establishes the main stages of health risk assessment associated with the migration of chemicals from materials intended for contact with food products: hazard identification, hazard characterization (allowing to determine the list of priority chemicals requiring further assessment, according to the proposed algorithm), exposure assessment (taking into account the determination of doses of substances entering the body with the diet), risk characterization; the calculation of approximate permissible amounts of migration of priority chemicals is given.

TECHNICAL ADVANTAGES

Technical advantages: prevention of harm to human health associated with the migration of chemical substances from packaging and materials intended for contact with food products.

Scientific and technical level: exceeds the level of the best domestic analogues.

EXPECTED RESULT OF APPLICATION

Practical use of the instructions for use in the practice of institutions carrying out state sanitary supervision and scientific health care organisations will contribute to improving the safety of materials through timely updating of both the lists of normative indicators and permissible levels of their migration; expanding the list of types of materials used in the food industry by assessing the health risks associated with the use of new types of materials; medical prevention of non-communicable diseases; and medical prevention of non-communicable diseases.

CURRENT STAGE OF DEVELOPMENT

At the implementation stage.

INFORMATION ON THE LEGAL PROTECTION

It was not planned.

POTENTIAL CONSUMERS AND INTERESTED PARTIES

Bodies and institutions exercising state sanitary supervision, and state medical research organizations authorized to conduct health risk assessments and establish state sanitary and hygienic standards for human environmental factors in the Republic of Belarus and abroad.

CONTACT INFORMATION

Vasily Tsygankov, Leading Researcher of the Laboratory of Systemic Problems of Food Hygiene, PhD, Associate Professor.

Phone number: (+375 17) 351 23 84

E-mail: pitanie_b@rspch.by

16. METHOD OF HYGIENIC STANDARDIZATION OF THE CONTENT OF ANTIBACTERIAL DRUGS IN THE AIR OF THE WORKING AREA

BRIEF DESCRIPTION OF INNOVATIVE DEVELOPMENT

Scope: the method can be used in a complex of medical services aimed at medical prevention of diseases and pathological conditions associated with the adverse effects of pharmaceutical substances of antibacterial drugs with bactericidal and/or bacteriostatic properties in the air of the working area.

Purpose: the implementation of the method will allow for the scientific substantiation of hygienic standards for antibacterial drugs in the air of the working area and will ensure a reduction in the potential risk of developing occupational allergic diseases to 1 % in pharmaceutical enterprises and healthcare institutions, and will allow for a reduction in the generalized indicator of the degree of allergy hazard to 4–5 points in pharmaceutical enterprises and healthcare institutions.

Main characteristics: method defines methodological approaches for the hygienic standardization of antibacterial drugs and ensures the rational development of maximum permissible concentrations of antibacterial drugs in the air of the working area, compliance with which during their production and use will ensure a reduction in the degree of allergy hazard in working conditions and the effective prevention of occupational allergies in workers in the pharmaceutical industry and healthcare.

TECHNICAL ADVANTAGES

Reduction of the degree of allergenic hazard in the working conditions of workers in the pharmaceutical industry and healthcare institutions from a potentially high level (up to 9 points in the absence of safe hygienic standards for antibiotics in the air of the working area) to a low-hazard level (4–5 points); reduction of the high potential risk of developing occupational allergic diseases in workers in the pharmaceutical industry and healthcare institutions (the probability of which in the absence of safe hygienic standards for antibiotics in the air of the working area can reach 30 %) to the minimum possible (probability up to 1 %).

Scientific and technical level: exceeds the level of the best domestic analogues.

EXPECTED RESULT OF APPLICATION

Practical use of the instructions in institutions implementing state sanitary supervision and scientific healthcare organizations will help reduce the degree of allergenic hazard in working conditions and reduce the potential risk of developing occupational allergic diseases among workers in the pharmaceutical industry and healthcare institutions.

CURRENT STAGE OF DEVELOPMENT

At the implementation stage.

INFORMATION ON THE LEGAL PROTECTION

It was not planned.

POTENTIAL CONSUMERS AND INTERESTED PARTIES

Healthcare organizations, including those implementing state sanitary supervision, state medical research organizations and other organizations implementing medical prevention, educational institutions with departments for the training, retraining and advanced training of specialists with higher education in the field of hygiene and preventive medicine and other organizations.

CONTACT INFORMATION

Ruslan Bogdanov, Head of the Industrial Toxicology Laboratory, Candidate of Medical Sciences, Associate Professor.

Phone number: (+ 375 17) 378 85 47

E-mail: promptox@rspch.by

17. METHODS OF HYGIENIC ASSESSMENT OF THE CHEMICAL COMPOSITION OF REPRESENTATIVE DIETS

BRIEF DESCRIPTION OF INNOVATIVE DEVELOPMENT

Scope of application: the methods can be used within a set of medical services aimed at medical prevention of diseases, as well as prevention of negative effects on health status associated with excessive or insufficient alimentary intake of chemical substances (nutrients and/or contaminants).

Purpose: the methods can be used for: identification of chemical substances entering through the alimentary route that are priority for hygienic assessment and hygienic regulation; assessment of public health risks associated with exposure to chemical substances in food products and diets (nutrients, chemical contaminants); substantiation of management decisions in the field of ensuring food safety and healthy nutrition of the population; development of programs for monitoring the safety of food products circulating on the market of the Republic of Belarus; assessment of the effectiveness of measures aimed at: reducing the level of contamination of food products with chemical contaminants and the associated alimentary exposure; optimization of qualitative and quantitative characteristics of nutrition of certain population groups.

Main characteristics: the methods make it possible to obtain objective data on the levels of chronic alimentary intake of chemical substances (nutrients and/or contaminants) as part of diets that generally correspond to the typical structure and consumption conditions.

TECHNICAL ADVANTAGES

Technical advantages: prevented damage to human health from the impact of the alimentary factor.

Scientific and technical level: corresponds to the levels of the best domestic and foreign analogues.

EXPECTED RESULT OF APPLICATION

The practical use of the methods in the practice of institutions exercising state sanitary supervision will contribute to improving the substantiation of management decisions in the field of ensuring food safety and healthy nutrition of the population, programs for laboratory monitoring of products circulating on the market of the Republic of Belarus, and the assessment of the effectiveness of measures aimed at reducing the level of contamination of food products with chemical contaminants and the associated alimentary exposure; it will contribute to the medical prevention of diseases associated with excessive or insufficient alimentary intake of chemical substances.

CURRENT STAGE OF DEVELOPMENT

Research or development (technological) work has been completed.

INFORMATION ON THE LEGAL PROTECTION

It was not planned.

POTENTIAL CONSUMERS AND INTERESTED PARTIES

Bodies and institutions exercising state sanitary supervision, and state medical research organizations authorized to conduct health risk assessments and establish state sanitary and hygienic standards for human environmental factors in the Republic of Belarus and abroad.

CONTACT INFORMATION

Ekaterina Fedorenko, Head of the International Activity Department, PhD, Associate Professor.

Phone number: (+ 375 17) 318 53 65

E-mail: fedorenko@rspch.by

IX. THE STATE INSTITUTION “REPUBLICAN SCIENTIFIC AND PRACTICAL CENTER FOR CHILDREN’S ONCOLOGY, HEMATOLOGY AND IMMUNOLOGY”

18. DEVELOPMENT OF A CAR RECEPTOR STRUCTURE FOR NK CELL MODIFICATION

BRIEF DESCRIPTION OF INNOVATIVE DEVELOPMENT

The development will be applied in the field of cellular immunotherapy of oncological diseases and is aimed at creating genetically modified NK cells. Its main purpose is to enhance the efficacy and specificity of the antitumor response of NK cells through (1) optimization of the structure of the CAR receptor itself (signaling domains) and (2) arming such a CAR receptor with cytokines (IL-15, 21, 18). These innovations are aimed at increasing the cytotoxic activity and persistence of NK cells in the tumor microenvironment.

Within the framework of the project, a library of next-generation chimeric antigen receptors (CARs) specifically adapted for NK cells is being created. Unlike standard CAR-T, which uses only CD3ζ, this development employs various combinations of domains: CD28, CD28TM, CD3z, 2B4, DAP12. The variability of their combinations will make it possible to determine the working activation and co-stimulatory signals and to avoid premature exhaustion of NK cells. In addition, arming CAR constructs with cytokines (IL-15, IL-21, IL-18 in secreted and membrane-bound forms) will achieve enhanced auto- and paracrine support of NK cells, also improving their proliferative activity and cytotoxicity without the systemic toxicity characteristic of administering recombinant cytokines.

TECHNICAL ADVANTAGES

Most clinical CAR-NK constructs use standard signaling domains (CD28, 4-1BB, CD3ζ) designed for T cells. The innovation of the project lies in the use of NK-specific domains — 2B4 and DAP12, — which more accurately reproduce the native signaling pathways of NK cells and reduce the risk of premature exhaustion. This distinguishes this development from those existing on the market. Similar developments are being carried out in the USA, China, South Korea, as well as Russia. However, publicly available data do not indicate the use of a combination of signaling domains with simultaneous arming with cytokines in a single lentiviral construct. This development is at the stage of creating a library of constructs and selecting the most suitable method for obtaining CAR-NK cells for the developer.

EXPECTED RESULT OF APPLICATION

It is expected to obtain CAR-NK cells (specific to the required antigen) with enhanced cytotoxicity, proliferative activity and resistance in the tumor microenvironment. This will make it possible to create the basis for a domestic targeted cell-based medicinal product, reduce dependence on foreign technologies and expand the possibilities of treating oncological diseases resistant to standard therapy.

CURRENT STAGE OF DEVELOPMENT

Research or development (technological) work has been completed.

INFORMATION ON THE LEGAL PROTECTION

It is planned to file an invention application with the National Center of Intellectual Property of the Republic of Belarus to protect the CAR receptor constructs and methods for modifying NK cells.

POTENTIAL CONSUMERS AND INTERESTED PARTIES

Institutions of the Ministry of Health, academic institutions, private enterprises of Belarus and other countries.

CONTACT INFORMATION

Yuliya Lebedeva, Junior Research Fellow, Postgraduate Student.

Phone number: (+375 29) 361 34 91

E-mail: lebedevaula859@gmail.com

19. BIOINFORMATICS ALGORITHM FOR AUTOMATED DETECTION AND INTERPRETATION OF GENETIC VARIANTS ASSOCIATED WITH PRIMARY IMMUNODEFICIENCIES

BRIEF DESCRIPTION OF INNOVATIVE DEVELOPMENT

The developed algorithm enables automated processing and interpretation of high-throughput sequencing data. The development is intended for the analysis of genetic variants identified during the examination of patients with suspected primary immunodeficiency.

The development is presented as a web resource with a graphical interface. It supports uploading sequencing data in FASTQ and VCF formats. Depending on the type of input data, a specific sequence of bioinformatics processing stages is performed, including quality assessment of sequencing data, mapping to a reference genome, detection and annotation of genetic variants, their comparison with information from specialized genetic databases, and generation of a final report in DOCX format.

The advantage of the development is the comprehensive automation of routine stages of bioinformatics processing and interpretation of high-throughput sequencing data, which reduces the amount of manual work for the specialist and shortens the time required to prepare the final conclusion. The software solution can be deployed locally and used with local copies of the necessary databases, ensuring operation under conditions of limited access to external internet resources.

The development is intended for use in research and medical organizations that employ next-generation sequencing technologies in the diagnosis of hereditary diseases.

TECHNICAL ADVANTAGES

Compared to foreign analogues, the development offers the following advantages: significantly lower cost of use and no need to purchase expensive licenses; ensuring the security of personal data through the possibility of local processing; reduced dependence on foreign software and information resources; and the ability to flexibly configure algorithms and analysis parameters according to the needs of a specific laboratory. No information has been identified about domestic analogues with comparable functionality.

EXPECTED RESULT OF APPLICATION

Implementation in the work of research and medical organizations of the country.

CURRENT STAGE OF DEVELOPMENT

A prototype sample has been produced.

INFORMATION ON THE LEGAL PROTECTION

Know-how.

POTENTIAL CONSUMERS AND INTERESTED PARTIES

Research and medical organizations using next-generation sequencing technologies in the diagnosis of hereditary diseases.

CONTACT INFORMATION

Irina Guryanova, Head of the Laboratory of Molecular Genetic Research of the Scientific Department.

Phone number: (+375 17) 287 10 69

E-mail: i.guryanova@oncology.by

Pavel Bobrik, Junior Researcher.

Phone number: (+375 33) 601 59 10

E-mail: bobrykpavel2001@gmail.com

X. STATE INSTITUTION “REPUBLICAN SCIENTIFIC AND PRACTICAL CENTER FOR PULMONOLOGY AND TUBERCULOSIS”

20. METHOD FOR DETERMINING THE DRUG SENSITIVITY OF *MYCOBACTERIUM ABSCESSUS* TO MACROLIDES AND AMINOGLYCOSIDES

BRIEF DESCRIPTION OF INNOVATIVE DEVELOPMENT

The method is based on the use of real-time polymerase chain reaction (PCR-RT) to determine the drug sensitivity of *Mycobacterium abscessus* to macrolides and aminoglycosides. The method allows for the rapid (within 6 hours) determination of the presence of drug sensitivity to macrolides and aminoglycosides, which is necessary when selecting a treatment regimen.

The method is implemented by sequentially completing the following stages: 1) preparation of biological material for research; 2) extraction of mycobacterial DNA from research material samples; 3) determination of macrolide resistance mutations in the *erm* and *rhl* genes and aminoglycoside resistance mutations in the *rrs* gene using PCR-RT; 4) recording and evaluation of the results. To implement the method, probes and primers specific to wild-type (FAM-labeled) and mutant alleles of the *erm(41)*, *rhl*, and *rrs* genes (JOE/R6G-labeled) have been developed. Fluorescence above the threshold level on the R6G/JOE and/or ROX channels indicates the presence of mutations responsible for the phenotypic resistance of *Mycobacterium abscessus* to macrolides and aminoglycosides.

The development is intended for laboratory diagnosticians, bacteriologists, and other specialist medical professionals in healthcare institutions who provide medical care to patients with diagnosed infections caused by *Mycobacterium* in inpatient and/or outpatient settings, including in day-care department settings. The field of application of the results obtained is phthisiology, pulmonology, and surgical diseases.

TECHNICAL ADVANTAGES

The method allows for the rapid (within 6 hours) determination of the presence of drug sensitivity to macrolides and aminoglycosides, which is essential when selecting a treatment regimen. The method has high sensitivity and specificity — 94 and 98 %, respectively — and makes it possible to identify samples containing a mixture of cultures. There are no domestic analogues. The molecular-genetic method used in the country to determine the drug sensitivity of non-tuberculous mycobacteria to macrolides and aminoglycosides is a commercial method that depends on imports.

EXPECTED RESULT OF APPLICATION

The economic advantages of the method lie in increasing the effectiveness of treating mycobacterial infections caused by *Mycobacterium abscessus* by selecting an appropriate treatment regimen for patients, in improving the quality of epidemiological studies on the prevalence of drug resistance among *Mycobacterium abscessus* complex mycobacteria, in strengthening the foundation of national diagnostic technologies, and in achieving import substitution.

CURRENT STAGE OF DEVELOPMENT

The instructions for use “Method for determining the drug sensitivity of *Mycobacterium abscessus* to macrolides and aminoglycosides”, registration No. 059-1125, approved 19.03.2026, have been developed and tested. The development is being implemented.

INFORMATION ON THE LEGAL PROTECTION

No.

POTENTIAL CONSUMERS AND INTERESTED PARTIES

Organizations providing medical care to patients with mycobacterial disease.

CONTACT INFORMATION

Liudmila Bohush, Leading Researcher at the Department of Pulmonology and Surgical Methods of Treating Diseases of the Respiratory System.
Phone number: (+375 29) 671 92 34
E-mail: ludabogush@mail.ru

21. ALGORITHM FOR SELECTING A CHEMOTHERAPY REGIMEN FOR PATIENTS WITH MULTIDRUG-RESISTANT TUBERCULOSIS, INCLUDING FLUOROQUINOLONE RESISTANCE

BRIEF DESCRIPTION OF INNOVATIVE DEVELOPMENT

The innovation is an algorithm for selecting and modifying chemotherapy regimens for patients with multidrug-resistant tuberculosis, including fluoroquinolone-resistant tuberculosis.

Its field of application is tuberculosis care.

The algorithm is intended for tuberculosis specialists and other physicians providing care to patients with drug-resistant tuberculosis.

Its purpose is to support a systematic, individualized selection of chemotherapy regimens based on the pathogen's drug susceptibility, previous treatment, comorbidities, contraindications, and drug tolerability.

The algorithm is structured as a multistage flowchart integrating baseline assessment, chemotherapy regimen selection, interim evaluation of effectiveness and safety, treatment adherence monitoring, and determination of subsequent treatment steps. It specifies the conditions for switching regimens in cases of treatment failure, detection of additional drug resistance, or adverse events.

Implementation relies on drug susceptibility testing using molecular genetic and culture-based methods, together with clinical, laboratory, electrocardiographic, imaging, and microbiological findings. These results and individual patient characteristics determine eligibility for a six-month, nine-month, or individualized longer chemotherapy regimen. For patients with HIV infection, the algorithm also considers the concurrent use of anti-tuberculosis and antiretroviral medicines.

Key operational characteristics include a clear visual structure, a consistent decision-making sequence, opportunities for reassessment and regimen modification throughout treatment, and applicability in inpatient and outpatient settings using healthcare facilities' laboratory and diagnostic resources.

The distinctive feature of the innovation is the integration of chemotherapy regimen selection and subsequent monitoring of treatment effectiveness, safety, and tolerability into a single algorithm.

TECHNICAL ADVANTAGES

The innovation integrates chemotherapy regimen selection and modification, taking into account the pathogen's drug susceptibility, treatment effectiveness, safety, tolerability, and the characteristics of the national patient cohort. It provides for six-month, nine-month, and individualized longer treatment regimens.

Treatment effectiveness with 24-week regimens was 91.2–94.6 %, comparable to outcomes reported in international clinical trials of similar short regimens (84–96 %). Sputum culture conversion by day 90 was achieved in 98.23 % of patients. Serious adverse events were reported in 11.5 % of patients.

The economic advantages of six-month regimens compared with individualized regimens lasting at least 18 months arise from shorter treatment duration, fewer required doses of medicines, and fewer outpatient visits. Modelled estimates indicate that wider use of shorter regimens could generate annual savings of up to USD 9.8 million and reduce the share of the national TB programme budget allocated to drug-resistant TB treatment from 25.8 to 12.6–12.8 %.

Additional socioeconomic benefits may include shorter periods of temporary work incapacity, lower expenditure on sickness benefits, and reduced economic losses associated with forgone GDP.

EXPECTED RESULT OF APPLICATION

Implementation of the developed algorithm is expected to support individualized selection and timely modification of chemotherapy regimens, improve treatment effectiveness and safety in patients with MDR-TB, including fluoroquinolone resistance, enhance treatment adherence, and reduce healthcare costs.

CURRENT STAGE OF DEVELOPMENT

Instructions for use of the chemotherapy regimen selection algorithm for patients with multidrug-resistant tuberculosis, including fluoroquinolone resistance, were approved by the Ministry of Health of the Republic of Belarus on October 30, 2025 (Registration No. 029-0625).

INFORMATION ON THE LEGAL PROTECTION

No.

POTENTIAL CONSUMERS AND INTERESTED PARTIES

Organizations providing TB care to patients with multidrug-resistant tuberculosis, including cases with fluoroquinolone resistance.

CONTACT INFORMATION

Natalia Yatskevich, Leading Researcher, Department of Laboratory Diagnostics and Treatment of Tuberculosis, RSPC for Pulmonology and Tuberculosis, MD, PhD, Associate Professor.

Phone number: (+375 29) 703 09 02

E-mail: yahoravanatallia@mail.ru

22. INSTRUCTION DESCRIBING THE “METHODS FOR IDENTIFICATION OF THE *MYCOBACTERIUM TUBERCULOSIS* GENOTYPE”

BRIEF DESCRIPTION OF INNOVATIVE DEVELOPMENT

The development is an instruction for the application of molecular genetic methods for identification of *Mycobacterium tuberculosis* genotypes not belonging to the Beijing genetic family. The method is intended for molecular epidemiological surveillance of tuberculosis, investigation of infection transmission, assessment of in-hospital spread of the pathogen, differentiation between relapse and reinfection, and evaluation of associations between genotype, drug resistance and clinical characteristics.

The method is based on real-time polymerase chain reaction (qPCR) using genotype-specific molecular markers and fluorescence-based detection of results. It enables identification of the major non-Beijing genotypes of *M. tuberculosis*, including LAM, Haarlem, T, Ural, X, S, Manu and EAI.

The technology can be implemented in laboratories equipped with standard real-time PCR instrumentation without the need for complex multistep genotyping procedures. Automated recording of results is possible, providing standardized and reproducible testing.

The development is intended for use in tuberculosis healthcare laboratories, research institutions, and units performing molecular epidemiological surveillance of tuberculosis.

TECHNICAL ADVANTAGES

Compared with conventional *M. tuberculosis* genotyping methods, including spoligotyping and MIRU-VNTR, the developed technology is based on real-time PCR and is characterized by lower labor intensity, fewer analytical steps, and reduced testing costs.

The technology provides high specificity and reproducibility of results, automated data recording, and standardized testing. The analysis can be performed using real-time PCR equipment already available in healthcare laboratories, without requiring specialized instrumentation.

The development expands existing molecular epidemiological surveillance, which is predominantly focused on identification of the Beijing genotype, by enabling identification of the major non-Beijing genotypes

of *M. tuberculosis*. This provides a more complete characterization of the circulating *M. tuberculosis* population and enables genotyping to be used directly in epidemiological investigations.

EXPECTED RESULT OF APPLICATION

Implementation of the development will enable identification of the major non-Beijing genotypes of *M. tuberculosis* at the genotype level, improve molecular epidemiological surveillance, facilitate timely investigation of infection transmission chains and in-hospital spread of the pathogen, and differentiate relapse from reinfection. The expected outcome is improved infection control and more objective characterization of the circulating *M. tuberculosis* population.

CURRENT STAGE OF DEVELOPMENT

Research or development (technological) work has been completed.

INFORMATION ON THE LEGAL PROTECTION

Patent protection is not available. The development has been formalized as the instruction "Methods for Identification of *Mycobacterium tuberculosis* Genotype", registration No. 006-1125.

POTENTIAL CONSUMERS AND INTERESTED PARTIES

In the Republic of Belarus: republican and regional tuberculosis control organizations, clinical diagnostic and microbiological laboratories, healthcare organizations performing molecular genetic testing, research institutions, and public health centers.

Abroad: tuberculosis control services and reference laboratories in countries with a high prevalence of tuberculosis, drug-resistant tuberculosis, and genetically diverse *M. tuberculosis* populations; research centers conducting molecular epidemiological surveillance of tuberculosis.

CONTACT INFORMATION

Veronica Slizen, PhD, Associate Professor.

Phone number: (+375 29) 260 32 68

E-mail: veronal75@yandex.ru

XI. REPUBLICAN RESEARCH AND PRACTICAL CENTER OF NEUROLOGY AND NEUROSURGERY, REPUBLICAN INNOVATIVE UNITARY ENTERPRISE “SCIENCE AND TECHNOLOGY PARK OF BNTU “POLYTECHNIC”

23. MEDICAL DEVICE “SYNTHETIC PROTECTOR FOR MICROVASCULAR DECOMPRESSION”

BRIEF DESCRIPTION OF INNOVATIVE DEVELOPMENT

The medical device (MD) “Synthetic Protector for Microvascular Decompression” has been created for use in surgical treatment — microvascular decompression of the facial, vestibulocochlear, and other cranial nerves. The developed Teflon protector, which has a specific shape and hydrophobicity, serves as a material for positional separation of the nerves and the adjacent vessel.

TECHNICAL ADVANTAGES

An import-substituting MD “Synthetic Protector for Microvascular Decompression” has been created; there are no analogues in the Republic of Belarus.

Intraoperatively, the protector is modeled in the form of an elongated soft roller or a set of small fragments that allow creating a multi-point support for the vascular loop. The thickness of the protector is selected to be minimal, sufficient only to separate the vessel and the nerve, without creating pressure on the structures of the facial-vestibular bundle.

The protector performs the role of a passive support that allows stabilizing the new position of the vessel. After mobilization, the vessel is moved anteriorly, superiorly, or posteriorly (depending on the type of conflict), and the protector ensures stable elimination of contact between the vessel and the VII–VIII nerves by resting against the dense elements of the surgical field and preventing its recurrent displacement toward the facial-vestibular bundle.

EXPECTED RESULT OF APPLICATION

The use of the MD “Synthetic Protector for Microvascular Decompression” will improve the quality of medical care for patients with hemi facial spasm and vestibular disorders caused by vascular compression of the vestibulocochlear nerve.

CURRENT STAGE OF DEVELOPMENT

The development has been introduced into production.

INFORMATION ON THE LEGAL PROTECTION

Registration certificate for “Synthetic Protector for Microvascular Decompression” TU BY 100232486.064-2025 No. IM-7.117479.

POTENTIAL CONSUMERS AND INTERESTED PARTIES

Neurosurgery departments of healthcare institutions of the Republic of Belarus and CIS countries.

CONTACT INFORMATION

Irina Maryenko, Chief Researcher of the Neurological Department of GI “Republican Research and Practical Center of Neurology and Neurosurgery”, Head of the Republican Center for Paroxysmal States, Doctor of Medical Sciences, Associate Professor.

Phone number: (+375 29) 653 10 50

E-mail: iramaryenko@gmail.com

Ryszard Sidorovich, Head of the Neurosurgical Department of GI “Republican Research and Practical Center of Neurology and Neurosurgery”, Doctor of Medical Sciences, Professor.

Phone number: (+375 29) 130 82 17

E-mail: info@neuro.by

Sample of a synthetic protector for microvascular decompression (in the form of a roll, 2 g)

XII. STATE INSTITUTION "REPUBLICAN SCIENTIFIC AND PRACTICAL CENTER FOR TRANSFUSIOLOGY AND MEDICAL BIOTECHNOLOGY"

24. THE DEVELOPMENT OF THE DOMESTIC MEDICINAL PRODUCT "PLAZMABEL" BASED ON VIRUS INACTIVATED CRYOPLASMA

Domestic medicinal product "PlasmaBel"

BRIEF DESCRIPTION OF INNOVATIVE DEVELOPMENT

A technology has been developed for producing the domestic medicinal product "PlasmaBel", which is solvent/detergent treated human plasma. Optimal conditions for producing the medicinal product have been selected at all stages of production, which ensures the preservation of a high level of factorial potential. Preclinical studies have confirmed the safety and non-toxicity of the medicinal product. Validation of the viral inactivation used in the production of the medicinal product confirms the correctness of the choice of the plasma treatment method with S/D reagents. The medicinal product has undergone full quality control in accordance with the State Pharmacopoeia of the Republic of Belarus. Clinical trials are ongoing.

TECHNICAL ADVANTAGES

There are no domestic analogues. It is a strategic medicinal product, the development of a production technology for which will ensure the national biosafety of the Republic of Belarus.

EXPECTED RESULT OF APPLICATION

Improving the quality of medical services.

CURRENT STAGE OF DEVELOPMENT

Research or development (technological) work is being carried out; pilot scale batches of the medicinal product have been produced, and their full control has been carried out in accordance with the regulatory documents.

INFORMATION ON THE LEGAL PROTECTION

No.

POTENTIAL CONSUMERS AND INTERESTED PARTIES

Neurosurgery departments of healthcare institutions of the Republic of Belarus and CIS countries.

CONTACT INFORMATION

Maryia Shcherbina, Researcher at the Laboratory of Industrial Transfusiology, Postgraduate Student.

Phone number: (+375 33) 677 00 81

E-mail: shcherbina@blood.by

XIII. THE EDUCATIONAL, SCIENTIFIC, AND PRODUCTION REPUBLICAN UNITARY ENTERPRISE "UNITECHPROM BSU", THE BSU INSTITUTION "RESEARCH INSTITUTE OF PHYSICO-CHEMICAL PROBLEMS"

25. ORIGINAL ANTITUMOR DRUG "PROSPIDELONG"

BRIEF DESCRIPTION OF INNOVATIVE DEVELOPMENT

The original drug "Prospidelong" is intended for the treatment of patients with gastric cancer with disseminated peritoneal lesions. "Prospidelong" is the world's first antitumor drug for intraperitoneal use, providing a long-term antitumor effect directly at the site of application. The drug is a powder of bioabsorbable gel-forming dextran phosphate with immobilized prospidin for the preparation of a gel for local chemotherapy.

TECHNICAL ADVANTAGES

The main advantages of the drug: targeted transport of the cytostatic into the tumor bed and the residual component of the removed tumor; prolongation of the therapeutic effect; potentiation of the effect of radiation therapy; ensuring homeostasis; complete resorption after implantation. In *in vitro* and *in vivo* experiments, the antiproliferative effect of the drug "Prospidelong" significantly exceeds the effect of prospidin.

EXPECTED RESULT OF APPLICATION

It is assumed that the intraoperative use of the drug "Prospidelong" as a means of long-term local chemotherapy in patients with gastric cancer with disseminated peritoneal lesions will increase the duration and improve the quality of life of patients, as well as increase the duration of the relapse-free period.

CURRENT STAGE OF DEVELOPMENT

Clinical trials are underway.

INFORMATION ON THE LEGAL PROTECTION

Trademarks for the drug prospidelong under application No. 20241254 dated 05/07/2024.

Patents of the Republic of Belarus: No. 14762, 15136.

Patents of the Russian Federation: No. 2455007, 2442586.

POTENTIAL CONSUMERS AND INTERESTED PARTIES

Clinical institutions of oncological profile.

CONTACT INFORMATION

Pavel Bychkovsky, Candidate of Chemical Sciences, Associate Professor.

Phone number: (+375 17) 209 58 40

E-mail: p.bychkovsky@zhivitsapharm.by


The original drug "Prospidelong"

XIV. EDUCATIONAL INSTITUTION “BELARUSIAN STATE MEDICAL UNIVERSITY”

26. SYNERGISTIC COMBINATION OF BEMITYL AND THYMOQUINONE: NEW HORIZONS IN THE FIGHT AGAINST HYPOXIA AND EXTREME CONDITIONS

BRIEF DESCRIPTION OF INNOVATIVE DEVELOPMENT

Hypoxic conditions caused by pulmonary edema (including the side effects of severe viral infections, such as COVID-19) or asphyxia under extreme conditions require immediate and highly effective pharmacological intervention. Existing monotherapies are often limited in efficacy or exhibit pronounced side effects when dosages are increased.

The effect of thymoquinone, bemetyl, and their combinations on the lifespan of mice under conditions of acute hypoxia with hypercapnia following intragastric administration ($M \pm m$)

Notes:

T — thymoquinone, B — bemetyl; doses are given in mg/kg;

* — $p < 0.05$ compared to placebo according to Duncan's criterion;

— $p < 0.05$ relative to the individual effect of the same dose of thymoquinone, Duncan's criterion.

We offer a patented synergistic combination of two components with fundamentally different yet complementary mechanisms of action:

1. Bemetyl is a benchmark domestic antihypoxant and actoprotector with proven efficacy, traditionally used in military medicine, disaster medicine, and high-performance sports.
2. Thymoquinone is a key bioactive component of black cumin oil, a dietary supplement widely recognized and approved worldwide for its antioxidant, anti-inflammatory, and immunomodulatory properties.

During preclinical studies on an experimental model of hypoxia in a hermetic volume (which closely simulates the conditions of suffocation, acute respiratory distress syndrome, and pulmonary edema), the combination demonstrated unique results:

- enhanced therapeutic effect: thymoquinone multiplies (potentiates) the antihypoxic activity and efficacy of bemetyl;
- mutual detoxification and safety: possessing a powerful detoxification potential, bemetyl neutralizes and minimizes any possible side effects of thymoquinone, making the combination as safe as possible for the body.

The development holds colossal potential for the pharmaceutical market and state security:

- 100 % Production localization: the synthesis of the pharmaceutical substance of “Bemetyl” has been fully mastered at the Institute of Physical Organic Chemistry of the National Academy of Sciences of Belarus;

- independence from imports: the production of black cumin oil (the source of thymoquinone) is currently being actively commercialized in our country; this guarantees an uninterrupted raw material base, cost accessibility, and high export potential.

Both components of the synergistic combination already possess official regulatory status, which dramatically accelerates and simplifies their integration into clinical practice:

- Military Medicine and EMERCOM: increasing the survival rate of personnel under conditions of oxygen deficiency, smoke inhalation, and toxic lesions;
- Public Health: inclusion in treatment protocols for pulmonological pathologies accompanied by hypoxia and pulmonary edema;
- Sports Medicine: legal (non-doping) enhancement of endurance and accelerated rehabilitation of athletes through the combined use of a proven actoprotector and a globally approved dietary supplement.

TECHNICAL ADVANTAGES

For the first time in global practice, the principle of the combined use of a synthetic actoprotector-antihypoxant with a plant-derived inducer of antioxidant defense has been implemented. The patent purity confirms the complete absence of direct foreign analogues with a similar mechanism of mutual enhancement (potentiation) of the effect.

The combination allows for a reduction in the effective dose of bemityl while maintaining (or increasing) its antihypoxic action, which drastically expands the therapeutic window and minimizes the toxic load on the liver and kidneys.

The use of ready-made bemityl, which has been fully mastered in production at the Institute of Physical Organic Chemistry of the National Academy of Sciences of Belarus, in combination with accessible plant raw materials, ensures that the cost of the final form is 3 to 5 times lower than that of foreign analogues of similar therapeutic focus.

EXPECTED RESULT OF APPLICATION

The implementation of the patented combination will provide:

- in medicine: a manifold increase in survival rates and body resilience to critical oxygen deficiency during pulmonary edema, asphyxia, and extreme conditions;
- in sports: legal acceleration of adaptation to ultimate physical workloads;
- in pharmacology: a reduction in effective doses of bemityl and minimization of toxic risks due to detoxification synergism.

The project guarantees 100 % import substitution and technological sovereignty through a full-cycle production of substances in the Republic of Belarus based on the Institute of Physical Organic Chemistry of the National Academy of Sciences of Belarus (bemityl) and farmeko.by (black cumin oil dietary supplement).

CURRENT STAGE OF DEVELOPMENT

Research or development (technological) work has been completed.

INFORMATION ON THE LEGAL PROTECTION

Patent for a synergistic combination BY 24695 / A. U. Vauchok, O. S. Rashkevich; patent issue date: 01.09.2025; publication date: 20.09.2025.

POTENTIAL CONSUMERS AND INTERESTED PARTIES

Healthcare organizations of the Republic of Belarus, sports organizations.

CONTACT INFORMATION

Aleksandr Vauchok, Deputy Head of the Department of Pharmacology, Candidate of Medical Sciences, Associate Professor.

Phone number: (+375 29) 667 31 42

E-mail: vauchok@bsmu.by

Olga Rashkevich, Senior Lecturer of the Department of Pharmacology, Master of Medical Sciences.


Phone number: (+375 29) 264 94 68

E-mail: ola.rashkevich@gmail.com

XV. EDUCATIONAL INSTITUTION “GOMEL STATE MEDICAL UNIVERSITY”

27. ANTIBACTERIAL TUBULAR VASCULAR IMPLANT

BRIEF DESCRIPTION OF INNOVATIVE DEVELOPMENT

Intended for use in vascular surgery during reconstructive procedures involving pathologically altered human blood vessels in the ileofemoral segment. The vascular implants are produced by modifying synthetic prostheses made of polyester filaments (woven, knitted, or straight) through the application of a polymer coating based on biologically active substances. Sterilization is achieved via gamma irradiation. The implants exhibit the following characteristics: antimicrobial properties (following impregnation in a 1 % vancomycin solution), zero surgical porosity, and enhanced coating biocompatibility.

TECHNICAL ADVANTAGES

The developed vascular implant with a biopolymer coating has the following advantages: lower cost, absence of toxicity, biocompatibility of the material, which is replaced by connective tissue, and a prolonged antibacterial effect (lasting at least 5 days, which is a critical period for prosthetic implantation infection).

Direct competitors in general are:

- textile vascular prostheses with an antibacterial coating based on gelatin; however, there is a known drawback of gelatin — possible allergic reactions and the rapid leaching of the antibiotic from the coating, leading to the loss of the protective antibacterial effect, as well as the high cost of this product;
- vascular prostheses with a silver coating, in addition to their high cost, provide a constant antibacterial effect, which may potentially cause cellular toxicity.

EXPECTED RESULT OF APPLICATION

The effectiveness of the development will lie in improving outcomes and minimizing the costs of treating patients with chronic obliterating diseases of the lower limb arteries.

CURRENT STAGE OF DEVELOPMENT

A prototype was released.

INFORMATION ON THE LEGAL PROTECTION

Patent.

POTENTIAL CONSUMERS AND INTERESTED PARTIES

Organizations of the Ministry of Health of the Republic of Belarus.

CONTACT INFORMATION

Yauheni Doroshko, Associate Professor at the Department of Surgical Diseases No. 1 with the Course of Cardiovascular Surgery, Candidate of Medical Sciences.

Phone number: (+375 44) 766 78 45


E-mail: ronaldy93@mail.ru

XVI. VITEBSK STATE ORDER OF PEOPLES' FRIENDSHIP MEDICAL UNIVERSITY

28. COMPREHENSIVE DIAGNOSIS OF METAL HYPERSENSITIVITY

BRIEF DESCRIPTION OF INNOVATIVE DEVELOPMENT

Comprehensive diagnosis of metal hypersensitivity involves a scoring assessment of the likelihood of metal allergy and an evaluation of the need for further provocative allergy testing. Provocative allergy testing includes skin patch tests using metal compounds found in medical devices (in a petrolatum base) and a sublingual-oral test with titanium oxide.

TECHNICAL ADVANTAGES

Comprehensive diagnostics for metal hypersensitivity in patients with acute and chronic intolerance reactions to metal-containing medical devices is being used for the first time in the Republic of Belarus and has no foreign counterparts registered in the country. The diagnostic accuracy of this comprehensive diagnostic method exceeds 80 %; sensitivity is above 70 %, and specificity is above 90 %.

EXPECTED RESULT OF APPLICATION

The widespread clinical application of the developed comprehensive diagnostic method for metal hypersensitivity will yield significant social and economic benefits by preventing adverse reactions and treatment complications, as well as by improving patients' quality of life.

CURRENT STAGE OF DEVELOPMENT

Two instructions for use have been registered:

– Method for performing skin patch testing with metal compounds / Aliakhnovich N. S., Ishchanka A. U., Kuntsevich Z. S., Khishova O. M., Shchurok I. N.; approved June 19, 2026; reg. No. 081-1225. — 7 pp.

– Method for performing a sublingual-oral test with titanium oxide / Alyakhnovich N. S., Yanchanka U. U., Ishchanka A. U., Goncharov A. E.; approved September 10, 2025; reg. No. 112-1124. — 8 pp.

The development is put into practical use at the Vitebsk Regional Clinical Hospital (Vitebsk) and in the educational process at the Vitebsk State Order of Peoples' Friendship Medical University (Vitebsk).

INFORMATION ON THE LEGAL PROTECTION

No.

POTENTIAL CONSUMERS AND INTERESTED PARTIES

Healthcare organizations of the Republic of Belarus providing medical care to patients suffering from allergic reactions to metal-containing medical devices.

CONTACT INFORMATION

Natalia Aliakhnovich, Associate Professor at the Department of Clinical Immunology and Allergology (including the Faculty of Advanced Training and Retraining), Candidate of Medical Sciences, Associate Professor.

Phone number: (+375 29) 212 51 24

E-mail: alyahnovich@bk.ru

29. DUAL-LUMEN NEPHROSTOMY CATHETER

BRIEF DESCRIPTION OF INNOVATIVE DEVELOPMENT

A dual-lumen nephrostomy catheter comprising a main body formed by a thin, flexible tube and a pigtail loop; the main body consists of an inflow lumen and an aspiration lumen featuring drainage openings located along the outer and inner contours of the distal curved section of the loop, as well as a rounded tip where the aspiration lumen terminates blindly, while the inflow lumen has an opening for a guidewire.

TECHNICAL ADVANTAGES

A distinctive feature of the proposed utility model is the catheter's dual-lumen design, with drainage openings located in the distal sections of the loop curve; this enables continuous-flow irrigation — using saline or antiseptics (in cases of extrarenal abscesses) combined with active aspiration of contents — as well as flushing of the catheter should one of its drainage channels become obstructed.

EXPECTED RESULT OF APPLICATION

The use of the proposed dual-lumen nephrostomy catheter enables active drainage and lavage of purulent processes within a dilated pelvicalyceal system caused by obstruction, as well as of purulent foci in the renal parenchyma and paranephric space associated with purulent-inflammatory kidney diseases; it also allows for monitoring the quality of the "irrigation-aspiration" process through visual inspection of the aspirated fluid and its volume. Cleansing of the purulent cavity is achieved more rapidly and effectively, which has a positive impact on the quality of treatment for the purulent-septic process.

CURRENT STAGE OF DEVELOPMENT


A prototype was released.

INFORMATION ON THE LEGAL PROTECTION

Patent for an invention.

POTENTIAL CONSUMERS AND INTERESTED PARTIES

Surgical and urological inpatient facilities.

Dual-lumen
nephrostomy catheter

CONTACT INFORMATION

Konstantin Kubrakov, Professor at the Department of Neurology and Neurosurgery,
 Doctor of Medical Sciences, Associate Professor.
 Phone number: (+375 29) 736 93 36
 E-mail: k-kubrakov@yandex.ru

30. REVERSE RETINACULOTOMY**BRIEF DESCRIPTION OF INNOVATIVE DEVELOPMENT**

The reverse retinaculotomy features a bow-shaped hook with rounded edges, measuring 12.5 by 8.5 mm externally and 3.5 by 3.5 mm internally, with a thickness of 1.5 mm. It has an S-shaped curve that bends toward the hook along the instrument's axis, and a four-sided handle with transverse grooves for grip.

TECHNICAL ADVANTAGES

A distinctive feature of the proposed utility model is the presence of a working part and an arrow-shaped hook with rounded edges, which facilitate atraumatic use of the instrument. The inner part of the model is sharpened at a 20-degree angle, allowing for isolated gripping and cutting of ligaments without damaging surrounding tissues. Additionally, it includes an S-shaped curve, bent toward the tool's axis in the direction of the hook, positioned between the handle and the working part to ensure comfortable alignment along the limb's axis.

EXPECTED RESULT OF APPLICATION

The use of this tool allows for an isolated subcutaneous dissection of the wrist ligament, preventing injury to surrounding structures, and can be performed without the need for intraoperative ultrasound or endoscopic navigation.

CURRENT STAGE OF DEVELOPMENT

Research or development (technological) work has been completed.

INFORMATION ON THE LEGAL PROTECTION


Utility Model Patent "Reverse Retinaculotomy" No. 13973 dated February 12, 2026.

POTENTIAL CONSUMERS AND INTERESTED PARTIES

Surgeons working in both inpatient and outpatient settings.

CONTACT INFORMATION

Konstantin Kubrakov, Professor at the Department of Neurology and Neurosurgery,
 Doctor of Medical Sciences, Associate Professor.
 Phone number: (+375 29) 736 93 36
 E-mail: k-kubrakov@yandex.ru

Reverse retinaculotomy

31. ANDROID MOBILE APP "LIVERSCORE"**BRIEF DESCRIPTION OF INNOVATIVE DEVELOPMENT**

In a physician's practice, when treating patients, they must calculate a large number of indices and indicators and complete a series of questionnaires to determine further treatment strategies. Calculating complex formulas for these indices is time-consuming, or online calculators are required.

Age	56	years
Albumin	45	g/L ↔
Bilirubin (initial)	234	μmol/L ↔
Bilirubin (day 7)	156	μmol/L ↔
Creatinine	100	μmol/L ↔
PT	15	sec

0.023 points
 Good prognosis. Scores <0.45 predict a 6-month survival of 85%.

Copy Results 
 Next Steps >>

The “LIVeScore” mobile app includes clinical scales and calculators incorporated into diagnostic and treatment protocols approved by the Ministry of Health of the Republic of Belarus. The app features the most frequently used calculations in gastroenterology, cardiology, and nephrology.

The uniqueness of the mobile app lies in its proprietary neural network, VRVrisk, based on the analysis of a patient’s clinical and laboratory parameters. The app also features a unique sarcopenia assessment module, which can be used for patients with both gastroenterology and other specialties (gerontology, endocrinology, etc.).

The infusion therapy and parenteral nutrition calculation module is based on guidelines used in physician practice. This module is a universal tool for physicians in anesthesiology and intensive care units, as well as other departments, and allows for the acceleration of routine calculations and improved quality of intensive care.

Furthermore, the app organizes current Ministry of Health of the Republic of Belarus guidelines for liver diseases by nosology.

TECHNICAL ADVANTAGES

The uniqueness of the mobile app lies in its proprietary neural network, VRVrisk, based on an analysis of the patient’s clinical and laboratory parameters. The neural network predicts the presence of grade 3 esophageal varices, including those complicated by bleeding, with 88.5 % accuracy and rules out the aforementioned conditions with 99.36 % accuracy, achieving an overall correct diagnosis rate of 96.88 %. The ROC curve area was 0.99. The developed neural network enables the detection of high-grade esophageal varices without endoscopic intervention with a high degree of certainty. Using the neural network, it will identify patients at high risk of bleeding and adjust subsequent patient management. There are no similar solutions available either in Belarus or abroad.

The app also features a unique sarcopenia assessment module, which can be used in patients with gastroenterological and other specialties (gerontological, endocrinological, etc.). The module includes a set of tests and a questionnaire to identify muscle mass and strength disorders. It has no equivalent among other mobile apps or online calculators.

The app is entirely based on the protocols of the Ministry of Health of the Republic of Belarus, which distinguishes it from similar apps localized to foreign protocols.

EXPECTED RESULT OF APPLICATION

Widespread implementation of the developed application will reduce the time spent on calculating routine patient indicators and improve the quality of treatment and diagnostic procedures.

CURRENT STAGE OF DEVELOPMENT

A working version of the mobile application has been released and is available for download.

INFORMATION ON THE LEGAL PROTECTION

It is planned to register the application in accordance with the established procedure.

POTENTIAL CONSUMERS AND INTERESTED PARTIES

The Ministry of Health of the Republic of Belarus, healthcare institutions (hospitals, polyclinics, hospitals, centers and departments of gastroenterology, cardiology, nephrology, anesthesiology and intensive care, gerontology and endocrinology departments), specialist doctors of relevant profiles, educational institutions and advanced training of medical workers.

CONTACT INFORMATION

Vyacheslav Pryshchepenko, Associate Professor of the Department of Propaedeutics of Internal Diseases, PhD, Associate Professor.

Phone number: (+375 29) 724 92 94

E-mail: prslava92@gmail.com

XVII. EDUCATIONAL INSTITUTION “BELARUSIAN STATE UNIVERSITY OF INFORMATICS AND RADIOELECTRONICS”, STATE INSTITUTION “REPUBLICAN SCIENTIFIC AND PRACTICAL CENTER FOR PULMONOLOGY AND PHTHIOLOGY”

32. SYSTEM FOR THE AUTOMATED DETERMINATION AND PERSONALIZED ADJUSTMENT OF OXYGEN THERAPY PARAMETERS IN REAL TIME BASED ON THE PATIENT'S CURRENT CONDITION

BRIEF DESCRIPTION OF INNOVATIVE DEVELOPMENT


The high prevalence of respiratory diseases leading to respiratory failure makes the task of personalizing and automating the selection and adjustment of oxygen therapy parameters particularly urgent, since oxygen therapy is the primary method for correcting hypoxemia and improving survival in patients with severe respiratory failure. However, current standard oxygen therapy protocols involve setting fixed parameters for the oxygen-enriched mixture — specifically flow rate and session duration. These parameters are typically adjusted manually, based on the experience and expertise of medical personnel. Yet clinical practice demonstrates that a patient's condition is dynamic; oxygenation needs fluctuate over time due to various factors, including physical activity, sleep, medication, environmental changes, and disease progression. A lack of timely adjustment can result in either ineffective treatment (due to persistent hypoxemia) or excessive oxygen delivery, which carries risks such as hyperoxia, toxic effects on lung tissue, and inefficient use of medical gas.

To enable real-time, personalized adjustment of oxygen therapy parameters based on the patient's current condition, a therapeutic and diagnostic biotechnical system has been developed. It comprises three modules:

- Diagnostic module: includes a set of sensors to monitor oxygen saturation (range: 50–100 %, accuracy ± 3 %), heart rate (range: 3–300 bpm, accuracy ± 3 %), and respiratory rate. The sensors are interfaced with a microcontroller that performs initial signal processing and digitization.

Personalized oxygen therapy system

An example of how a personalized oxygen therapy system works

– Control and communication module: implements the adaptive adjustment algorithm. It comprises a processor, memory, and input/output interfaces, as well as a wireless communication unit for data transmission to a remote server (telemedicine channel).

– Actuator module: contains an oxygen source (concentrator, cylinder, or centralized supply line), adjustable flow-control valves (1–5 L/min), and a pneumatic mixing and humidification circuit.

Operating principle of the proposed system: the diagnostic module continuously measures the patient's physiological parameters and transmits them to the control and communication module; the latter compares the received values against target values and, in accordance with the developed methodology, generates commands for the actuator module. The actuator module then adjusts the control elements to deliver the oxygen mixture with the required parameters. This closes the biological feedback loop, enabling control based on the patient's condition.

TECHNICAL ADVANTAGES

The modular architecture, based on the principles of complex system design, ensures flexibility and scalability. The system can operate in various configurations: as a fully automated inpatient system, as a home-use device with remote monitoring capabilities, and as a portable monitor for screening respiratory diseases and sleep-disordered breathing.

The novelty lies in the automated, personalized selection of oxygen therapy parameters and the real-time adjustment of these values via a feedback loop. There are no direct analogues. Currently, a physician performs the selection and adjustment of therapy settings manually.

Compliance with industry standards:

- Ministry of Health of the Republic of Belarus Registration Certificate No. IM-7.108117/2506 dated June 5, 2025;
- Certificate of Approval of Measuring Instrument Type No. 18518, issued by the State Committee for Standardization of the Republic of Belarus on March 10, 2025;
- ISO 13485:2016 Quality Management System Certificate of Conformity.

Complies with safety standards:

- STB IEC 60601-1-2-2006;
- STB EN 55011-2006;
- GOST 30324.0-95.

EXPECTED RESULT OF APPLICATION

Social impact:

- early detection of respiratory diseases complicated by respiratory failure (COPD, sleep apnea, alveolitis, coronavirus infection, bronchial asthma, etc.);
- improved treatment efficacy and safety: fewer side effects, less disability and mortality, longer life;
- increased efficacy and safety of home-based therapy.

Economic impact:

- reduced burden on and costs to healthcare and social protection systems associated with frequent and prolonged hospitalizations and early disability in patients;
- reduced consumption of oxygen (by at least 60 %, $p < 0.001$), medications, and electricity;
- savings in the medical personnel working time;
- reduced length of patients' hospital stay.

CURRENT STAGE OF DEVELOPMENT

Development is introduced into production.

INFORMATION ON THE LEGAL PROTECTION

Patents of the Republic of Belarus for a utility model, industrial design and trademark have been received.

POTENTIAL CONSUMERS AND INTERESTED PARTIES

Areas of application:

- healthcare facilities;
- sanatoriums and health resorts;
- patients with respiratory failure (due to COVID-19, COPD, bronchial asthma, alveolitis, ALS, etc.) during the course of the disease and rehabilitation — both in inpatient settings and at home.

CONTACT INFORMATION

Oleg Zelmansky, Associate Professor of the Department of Information Security, PhD, Associate Professor.

Phone number: (+375 29) 390 09 14

E-mail: 7650772@rambler.ru

XVIII. APEX-REAL LIMITED LIABILITY COMPANY (RESIDENT OF BREST SCIENCE AND TECHNOLOGY PARK)

33. DEVELOPMENT AND IMPLEMENTATION OF A NEW ELECTRONIC SYSTEM FOR MANAGING AND DISTRIBUTING PATIENT FLOWS

BRIEF DESCRIPTION OF INNOVATIVE DEVELOPMENT

The objective of the project is the development and implementation of a new electronic system for managing and distributing patient flows in transfusion centers (departments), blood transfusion stations, and healthcare institutions of the Republic of Belarus.

The goal of this project is to create an effective electronic system for distributing and managing patient flows that will minimize patient waiting time and improve the organization of their flows.

The system will cover the full service cycle and support the management of several queues simultaneously (registration, consultation, plasma, blood, platelets, apheresis, procedure room, information desk) with configurable priority and activity for each queue.

TECHNICAL ADVANTAGES

Existing solutions, including both analogue and modern electronic queue management systems, often do not take into account the specific operational characteristics of particular transfusion institutions. For the most part, such solutions do not address the tasks assigned to effective queue management and the distribution of people flows, which creates the need for developing individual approaches.

The proposed system is designed from the outset around the actual donor routes: the ticket contains a step-by-step route indicating floors and rooms (registration, blood test, appointment with a specialist doctor, buffet, barcode scanning, donation of blood or its components, payment of monetary compensation), which allows the patient to navigate the institution independently without staff assistance.

EXPECTED RESULT OF APPLICATION

The implementation of the system makes it possible to minimize patient waiting time, streamline their flows and increase the throughput capacity of the institution by eliminating intersections between different categories of donors (blood, plasma, platelets, apheresis). The donor receives a clear step-by-step route, while the staff receives an up-to-date picture of room and queue loads, which reduces the burden on the registration desk and doctors. According to market monitoring data for the Republic of Belarus, the system has no analogues in transfusion centers (departments). The preliminary estimate of the market volume is approximately 250–300 thousand Belarusian rubles per year.

CURRENT STAGE OF DEVELOPMENT

Research and development (engineering/technological) work has been completed.

INFORMATION ON THE LEGAL PROTECTION

The software of the electronic queue management system (EQMS) must be registered in accordance with the established procedure with the National Center of Intellectual Property of the Republic of Belarus and must have a license certificate for use, issued by the rights holder and transferred to customers.

POTENTIAL CONSUMERS AND INTERESTED PARTIES

Regional transfusion centers; district and city transfusion departments; blood transfusion stations of the Republic of Belarus.

CONTACT INFORMATION

Andrey Sytko, Director.

Phone number: (+375 29) 726 91 98

E-mail: apeksreal@gmail.com

XIX. FEDERAL STATE BUDGETARY EDUCATIONAL INSTITUTION OF HIGHER EDUCATION “KEMEROVO STATE MEDICAL UNIVERSITY” OF THE MINISTRY OF HEALTH OF THE RUSSIAN FEDERATION

34. A DEVICE FOR SELF-COLLECTION OF VAGINAL DISCHARGE — THE “JUNO” PROBE

BRIEF DESCRIPTION OF INNOVATIVE DEVELOPMENT

The device for self-collection of vaginal discharge is a plastic rod with a tip and a holder. The plastic rod is equipped with a stopper in the form of semi-oval “wings”; the tip is made of a soft, easily absorbent material; the holder is in the form of a flat ring. The device is intended for self-collection of vaginal discharge and subsequent testing of the obtained material for high-oncogenic-risk HPV. The device is in individual sterile packaging with a container for the tip.

The field of application is gynecology (within the framework of cervical cancer screening, specifically HPV testing). Almost 100 % of cervical cancer is associated with high-oncogenic-risk HPV infection. Timely identification and classification of such women into a risk group will significantly improve both the number of screening participants and its quality.

TECHNICAL ADVANTAGES

The advantage of this development is the absence of domestic analogues. There is only one imported device on the market. Due to the design features of the device and the presence of an industrial partner, the cost will be two times lower than that of the imported analogue. Implementation in screening programs in the regions will make it possible to use the device for both urban and rural populations as part of HPV testing within the secondary prevention program for cervical cancer.

EXPECTED RESULT OF APPLICATION

Thanks to the use of a device for self-collection of vaginal discharge, it is possible to expand the coverage of the cervical cancer screening program for the female population. The results of HPV testing for high-oncogenic genotypes after a doctor's collection of the material and self-collection of the material are comparable, which has been confirmed by numerous studies. According to clinical guidelines regarding the secondary prevention of cervical cancer, HPV testing should be carried out for all women aged 30–65. As of January 1, 2026, there are more than 43 million women in this age group, of whom about 9 million live in rural areas. Using this technology, it is possible to reduce the incidence and mortality from cervical cancer.

CURRENT STAGE OF DEVELOPMENT

A prototype was released.

INFORMATION ON THE LEGAL PROTECTION

Utility model patent RU 222041 U1, 07.12.2023. Application No. 2023124506 dated 22.09.2023.

Utility model patent RU 222091 U1, 11.12.2023. Application No. 2023124507 dated 22.09.2023.

Patents No. 13600 and No. 13601 have been obtained in the Republic of Belarus

POTENTIAL CONSUMERS AND INTERESTED PARTIES

Sales via marketplaces, pharmacies, and the supply of test systems to regional healthcare facilities for HPV testing as part of a cervical cancer screening program, including for women in remote areas.

CONTACT INFORMATION

Natalya Artymuk, Head of the Department of Obstetrics and Gynecology
named after Professor G. A. Ushakova, Doctor of Medical Sciences, Professor.

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© ГКНТ, 2026